

ORTOPEDSKA BOLNIŠNICA VALDOLTRA
Jadranska cesta 31
6280 Ankaran
Slovenia

Declaration Letter MRI Safety Evaluation for TENDER “Sklepne endoproteze in implantati za kirurgijo spinalnih deformacij”

To whom it may concern,

The information provided below is only applicable to the devices listed in Appendix 1, which are part of the following brands:

- Corail®
- Pinnacle®
- Trilock®
- REEF™
- Sigma® Total Knee System
- Sigma® HP Partial Knee
- Attune® Total Knee System

Based on internal MRI test data, the product codes listed in Appendix 1, can be considered as MRI Conditional in a 3.0T or lower power MR scanner, per ASTM F2503-13.

Please note that DePuy Synthes’ MRI Safety Evaluation project is still ongoing. Per DePuy Synthes agreement with BSI, an MRI Safety Brochure has been made available to Healthcare Professionals. The results of the MRI Safety Evaluation are published in the MRI Safety Brochure, and as new data becomes available, the MRI Safety Brochure is updated. Please note that DePuy Synthes is currently in the process of updating the current version of the Brochure, to include the most recent results, including some the product codes listed in Appendix 1.

Updates to DePuy Synthes Labelling (Instructions for Use and Packaging information) to include specific wording from the MRI testing and ASTM symbols will be made following the appropriate regulatory submissions.

In the meantime, individual requests regarding MRI safety of specific knee or hip constructs can be made through the Medical Information Request (MIR) process of our Medical Affairs department. Per the procedure, the requests can be made by Healthcare Professionals regarding information for the on-label or off-label use of all DePuy Synthes and other J&J Medical Devices products. The requests can be made online at www.jjmedir.com.

Signed on behalf of DePuy Orthopaedics, Inc.

Signed: _____

Nicola Keedy

Regulatory Affairs Associate Director

Date: see electronic signature

Appendix 1 – Scope of codes

Product code	Product description
3L92507	Corail STD Size 8
3L92509	Corail STD Size 9
3L92510	Corail STD Size 10
3L92511	Corail STD Size 11
3L92512	Corail STD Size 12
3L92513	Corail STD Size 13
3L92514	Corail STD Size 14
3L92515	Corail STD Size 15
3L92516	Corail STD Size 16
3L92518	Corail STD Size 18
3L92520	Corail STD Size 20
3L92498	Corail STD Collared Size 8
3L92499	Corail STD Collared Size 9
3L92500	Corail STD Collared Size 10
3L92501	Corail STD Collared Size 11
3L92502	Corail STD Collared Size 12
3L92503	Corail STD Collared Size 13
3L92504	Corail STD Collared Size 14
3L92505	Corail STD Collared Size 15
3L92506	Corail STD Collared Size 16
3L92508	Corail STD Collared Size 18
3L92521	Corail STD Collared Size 20
L20309	Corail NON Collared HO Size 9
L20310	Corail NON Collared HO Size 10
L20311	Corail NON Collared HO Size 11
L20312	Corail NON Collared HO Size 12
L20313	Corail NON Collared HO Size 13
L20314	Corail NON Collared HO Size 14
L20315	Corail NON Collared HO Size 15
L20316	Corail NON Collared HO Size 16
L20318	Corail NON Collared HO Size 18
L20320	Corail NON Collared HO Size 20
L971109	Corail Collared HO Size 9
L971110	Corail Collared HO Size 10
L971111	Corail Collared HO Size 11
L971112	Corail Collared HO Size 12
L971113	Corail Collared HO Size 13
L971114	Corail Collared HO Size 14

L971115	Corail Collared HO Size 15
L971116	Corail Collared HO Size 16
L971118	Corail Collared HO Size 18
L971120	Corail Collared HO Size 20
3L93709	Corail LAT COXA VARA Size 9
3L93710	Corail LAT COXA VARA Size 10
3L93711	Corail LAT COXA VARA Size 11
3L93712	Corail LAT COXA VARA Size 12
3L93713	Corail LAT COXA VARA Size 13
3L93714	Corail LAT COXA VARA Size 14
3L93715	Corail LAT COXA VARA Size 15
3L93716	Corail LAT COXA VARA Size 16
3L93718	Corail LAT COXA VARA Size 18
3L93720	Corail LAT COXA VARA Size 20
L971308	Corail Short Neck Collared size 8
L971309	Corail Short Neck Collared size 9
L971310	Corail Short Neck Collared size 10
L981308	Corail Short Neck NON Collared size 8
L981309	Corail Short Neck NON Collared size 9
L981310	Corail Short Neck NON Collared size 10
L971208	Corail Short Neck Collared size 8
L971209	Corail Short Neck Collared size 9
L971210	Corail Short Neck Collared size 10
L981208	Corail Short Neck NON Collared size 8
L971209	Corail Short Neck NON Collared size 9
L971210	Corail Short Neck NON Collared size 10
L20106	TIGE CORAIL2 DYSP S/A T6
L20006	TIGE CORAIL2 DYSP A/A T6
L98010	Corail Revision stem STD 10
L98011	Corail Revision stem STD 11
L98012	Corail Revision stem STD 12
L98013	Corail Revision stem STD 13
L98014	Corail Revision stem STD 14
L98015	Corail Revision stem STD 15
L98016	Corail Revision stem STD 16
L98018	Corail Revision stem STD 18
L98020	Corail Revision stem STD 20
L98110	Corail Revision stem HO 10
L98111	Corail Revision stem HO 11
L98112	Corail Revision stem HO 12
L98113	Corail Revision stem HO 13
L98114	Corail Revision stem HO 14

L98115	Corail Revision stem HO 15
L98116	Corail Revision stem HO 16
L98118	Corail Revision stem HO 18
L98120	Corail Revision stem HO 20
121701048	Pinnacle™ 100 Acetabular Cup 48mm
121701050	Pinnacle™ 100 Acetabular Cup 50mm
121701052	Pinnacle™ 100 Acetabular Cup 52mm
121701054	Pinnacle™ 100 Acetabular Cup 54mm
121701056	Pinnacle™ 100 Acetabular Cup 56mm
121701058	Pinnacle™ 100 Acetabular Cup 58mm
121701060	Pinnacle™ 100 Acetabular Cup 60mm
121701062	Pinnacle™ 100 Acetabular Cup 62mm
121701064	Pinnacle™ 100 Acetabular Cup 64mm
121701066	Pinnacle™ 100 Acetabular Cup 66mm
121720038	Pinnacle™ Sector Acetabular Cup 38mm
121720040	Pinnacle™ Sector Acetabular Cup 40mm
121720042	Pinnacle™ Sector Acetabular Cup 42mm
121720044	Pinnacle™ Sector Acetabular Cup 44mm
121720046	Pinnacle™ Sector Acetabular Cup 46mm
121722048	Pinnacle™ Sector Acetabular Cup 48mm
121722050	Pinnacle™ Sector Acetabular Cup 50mm
121722052	Pinnacle™ Sector Acetabular Cup 52mm
121722054	Pinnacle™ Sector Acetabular Cup 54mm
121722056	Pinnacle™ Sector Acetabular Cup 56mm
121722058	Pinnacle™ Sector Acetabular Cup 58mm
121722060	Pinnacle™ Sector Acetabular Cup 60mm
121722062	Pinnacle™ Sector Acetabular Cup 62mm
121722064	Pinnacle™ Sector Acetabular Cup 64mm
121722066	Pinnacle™ Sector Acetabular Cup 66mm
121922038	22,225mm Neutral Marathon™ Polyethylene liner OD 38mm
121922040	22,225mm Neutral Marathon™ Polyethylene liner OD 40mm
121922042	22,225mm Neutral Marathon™ Polyethylene liner OD 42mm
121928044	28mm Neutral Marathon™ Polyethylene liner OD 44mm
121928046	28mm Neutral Marathon™ Polyethylene liner OD 46mm
121928048	28mm Neutral Marathon™ Polyethylene liner OD 48mm
121928050	28mm Neutral Marathon™ Polyethylene liner OD 50mm
121928052	28mm Neutral Marathon™ Polyethylene liner OD 52mm
121928054	28mm Neutral Marathon™ Polyethylene liner OD 54mm
121928056	28mm Neutral Marathon™ Polyethylene liner OD 56mm

121928058	28mm Neutral Marathon™ Polyethylene liner OD 58mm
121928060	28mm Neutral Marathon™ Polyethylene liner OD 60mm
121928062	28mm Neutral Marathon™ Polyethylene liner OD 62mm
121928064	28mm Neutral Marathon™ Polyethylene liner OD 64mm
121928066	28mm Neutral Marathon™ Polyethylene liner OD 66mm
121932052	32mm Neutral Marathon™ Polyethylene Liner OD 52mm
121932054	32mm Neutral Marathon™ Polyethylene Liner OD 54mm
121932056	32mm Neutral Marathon™ Polyethylene Liner OD 56mm
121932058	32mm Neutral Marathon™ Polyethylene Liner OD 58mm
121932060	32mm Neutral Marathon™ Polyethylene Liner OD 60mm
121932062	32mm Neutral Marathon™ Polyethylene Liner OD 62mm
121932064	32mm Neutral Marathon™ Polyethylene Liner OD 64mm
121932066	32mm Neutral Marathon™ Polyethylene Liner OD 66mm
121936056	36mm Neutral Marathon™ Polyethylene Liner OD 56mm
121936058	36mm Neutral Marathon™ Polyethylene Liner OD 58mm
121936060	36mm Neutral Marathon™ Polyethylene Liner OD 60mm
121936062	36mm Neutral Marathon™ Polyethylene Liner OD 62mm
121936064	36mm Neutral Marathon™ Polyethylene Liner OD 64mm
121936066	36mm Neutral Marathon™ Polyethylene Liner OD 66mm
121883744	28mm Ceramax™ Alternative Bearing Insert OD 44mm
121883746	28mm Ceramax™ Alternative Bearing Insert OD 46mm
121883748	28mm Ceramax™ Alternative Bearing Insert OD 48mm
121883750	28mm Ceramax™ Alternative Bearing Insert OD 50mm
121883752	28mm Ceramax™ Alternative Bearing Insert OD 52mm
121883754	28mm Ceramax™ Alternative Bearing Insert OD 54mm
121883756	28mm Ceramax™ Alternative Bearing Insert OD 56mm
121883758	28mm Ceramax™ Alternative Bearing Insert OD 58mm
121883760	28mm Ceramax™ Alternative Bearing Insert OD 60mm
121883762	28mm Ceramax™ Alternative Bearing Insert OD 62mm
121883764	28mm Ceramax™ Alternative Bearing Insert OD 64mm
121883766	28mm Ceramax™ Alternative Bearing Insert OD 66mm
121882748	32mm Ceramax™ Alternative Bearing Insert OD 48mm
121882750	32mm Ceramax™ Alternative Bearing Insert OD 50mm

121882752	32mm Ceramax™ Alternative Bearing Insert OD 52mm
121881752	36mm Ceramax™ Alternative Bearing Insert OD 52mm
121881754	36mm Ceramax™ Alternative Bearing Insert OD 54mm
121881756	36mm Ceramax™ Alternative Bearing Insert OD 56mm
121881758	36mm Ceramax™ Alternative Bearing Insert OD 58mm
121881760	36mm Ceramax™ Alternative Bearing Insert OD 60mm
121881762	36mm Ceramax™ Alternative Bearing Insert OD 62mm
121881764	36mm Ceramax™ Alternative Bearing Insert OD 64mm
121881766	36mm Ceramax™ Alternative Bearing Insert OD 66mm
121928148	Marathon liner +4 10D 28mm/48
121928150	Marathon liner +4 10D 28mm/50
121928152	Marathon liner +4 10D 28mm/52
121928154	Marathon liner +4 10D 28mm/54
121928156	Marathon liner +4 10D 28mm/56
121928158	Marathon liner +4 10D 28mm/58
121928160	Marathon liner +4 10D 28mm/60
121928162	Marathon liner +4 10D 28mm/62
121928164	Marathon liner +4 10D 28mm/64
121928166	Marathon liner +4 10D 28mm/66
121932148	Marathon liner +4 10D 32mm/48
121932150	Marathon liner +4 10D 32mm/50
121932152	Marathon liner +4 10D 32mm/52
121932154	Marathon liner +4 10D 32mm/54
121932156	Marathon liner +4 10D 32mm/56
121932158	Marathon liner +4 10D 32mm/58
121932160	Marathon liner +4 10D 32mm/60
121932162	Marathon liner +4 10D 32mm/62
121932164	Marathon liner +4 10D 32mm/64
121932166	Marathon liner +4 10D 32mm/66
121936152	Marathon liner +4 10D 36mm/52
121936154	Marathon liner +4 10D 36mm/54
121936156	Marathon liner +4 10D 36mm/56
121936158	Marathon liner +4 10D 36mm/58
121936160	Marathon liner +4 10D 36mm/60
121936162	Marathon liner +4 10D 36mm/62
121936164	Marathon liner +4 10D 36mm/64
121936166	Marathon liner +4 10D 36mm/66
136529000	22,225mm 12/14 Articul/eze Ultamet™ Head 4
136530000	22,225mm 12/14 Articul/eze Ultamet™ Head 7
136511000	28mm 12/14 Articul/eze Ultamet™ Head 1.5 GR

136512000	28mm 12/14 Articul/eze Ultamet™ Head 5 BR
136513000	28mm 12/14 Articul/eze Ultamet™ Head 8.5 BL
136514000	28mm 12/14 Articul/eze Ultamet™ Head 12 BLK
136521000	32mm 12/14 Articul/eze Ultamet™ Head 1
136522000	32mm 12/14 Articul/eze Ultamet™ Head 5
136523000	32mm 12/14 Articul/eze Ultamet™ Head 9
136524000	32mm 12/14 Articul/eze Ultamet™ Head 13
136550000	36mm 12/14 Articul/eze Ultamet™ Head -2
136551000	36mm 12/14 Articul/eze Ultamet™ Head 1,5
136552000	36mm 12/14 Articul/eze Ultamet™ Head 5
136553000	36mm 12/14 Articul/eze Ultamet™ Head 8.5
136554000	36mm 12/14 Articul/eze Ultamet™ Head 12
136528310	28mm 12/14 Articul/eze Ceramax™ Head 1,5
136528320	28mm 12/14 Articul/eze Ceramax™ Head 5
136528330	28mm 12/14 Articul/eze Ceramax™ Head 8,5
136532310	32mm 12/14 Articul/eze Ceramax™ Head 1
136532320	32mm 12/14 Articul/eze Ceramax™ Head 5
136532330	32mm 12/14 Articul/eze Ceramax™ Head 9
136536310	36mm 12/14 Articul/eze Ceramax™ Head 1,5
136536320	36mm 12/14 Articul/eze Ceramax™ Head 5
136536330	36mm 12/14 Articul/eze Ceramax™ Head 8.5
136536340	36mm 12/14 Articul/eze Ceramax™ Head 12
121715500	PINN.6.5 CANCELLOUS SCREW 15mm
121720500	PINN.6.5 CANCELLOUS SCREW 20mm
121725500	PINN.6.5 CANCELLOUS SCREW 25mm
121730500	PINN.6.5 CANCELLOUS SCREW 30mm
121735500	PINN.6.5 CANCELLOUS SCREW 35mm
121740500	PINN.6.5 CANCELLOUS SCREW 40mm
121745500	PINN.6.5 CANCELLOUS SCREW 45mm
121750500	PINN.6.5 CANCELLOUS SCREW 50mm
121755500	PINN.6.5 CANCELLOUS SCREW 55mm
121760500	PINN.6.5 CANCELLOUS SCREW 60mm
121765500	PINN.6.5 CANCELLOUS SCREW 65mm
121770500	PINN.6.5 CANCELLOUS SCREW 70mm
124603000	Pinnacle® Apex Hole Eliminator
101204005	Trilock® STD Offset size 0
101204010	Trilock® STD Offset size 1
101204020	Trilock® STD Offset size 2
101204030	Trilock® STD Offset size 3
101204040	Trilock® STD Offset size 4
101204050	Trilock® STD Offset size 5
101204060	Trilock® STD Offset size 6
101204070	Trilock® STD Offset size 7

101204080	Trilock® STD Offset size 8
101204090	Trilock® STD Offset size 9
101204100	Trilock® STD Offset size 10
101204110	Trilock® STD Offset size 11
101204120	Trilock® STD Offset size 12
101214005	Trilock® LAT. size 0
101214010	Trilock® LAT. size 1
101214020	Trilock® LAT. size 2
101214030	Trilock® LAT. size 3
101214040	Trilock® LAT. size 4
101214050	Trilock® LAT. size 5
101214060	Trilock® LAT. size 6
101214070	Trilock® LAT. size 7
101214080	Trilock® LAT. size 8
101214090	Trilock® LAT. size 9
101214100	Trilock® LAT. size 10
101214110	Trilock® LAT. size 11
101214120	Trilock® LAT. size 12
121730038	Pinnacle Sector Gription 38mm
121730040	Pinnacle Sector Gription 40mm
121730042	Pinnacle Sector Gription 42mm
121730044	Pinnacle Sector Gription 44mm
121730046	Pinnacle Sector Gription 46mm
121732048	Pinnacle Sector Gription 48mm
121732050	Pinnacle Sector Gription 50mm
121732052	Pinnacle Sector Gription 52mm
121732054	Pinnacle Sector Gription 54mm
121732056	Pinnacle Sector Gription 56mm
121732058	Pinnacle Sector Gription 58mm
121732060	Pinnacle Sector Gription 60mm
121732062	Pinnacle Sector Gription 62mm
121732064	Pinnacle Sector Gription 64mm
121732066	Pinnacle Sector Gription 66mm
122128044	ALTRX Neutral 28x44mm
122128046	ALTRX Neutral 28x46mm
122128048	ALTRX Neutral 28x48mm
122128050	ALTRX Neutral 28x50mm
122128052	ALTRX Neutral 28x52mm
122128054	ALTRX Neutral 28x54mm
122132048	ALTRX Neutral 32x48mm
122132050	ALTRX Neutral 32x50mm
122132052	ALTRX Neutral 32x52mm
122132054	ALTRX Neutral 32x54mm

122132056	ALTRX Neutral 32x56mm
122132058	ALTRX Neutral 32x58mm
122136052	ALTRX Neutral 36x52mm
122136054	ALTRX Neutral 36x54mm
122136056	ALTRX Neutral 36x56mm
122136058	ALTRX Neutral 36x58mm
122136060	ALTRX Neutral 36x60mm
122136062	ALTRX Neutral 36x62mm
122136064	ALTRX Neutral 36x64mm
122136066	ALTRX Neutral 36x66mm
101204005	Trilock® STD Offset size 0
101204010	Trilock® STD Offset size 1
101204020	Trilock® STD Offset size 2
101204030	Trilock® STD Offset size 3
101204040	Trilock® STD Offset size 4
101204050	Trilock® STD Offset size 5
101204060	Trilock® STD Offset size 6
101204070	Trilock® STD Offset size 7
101204080	Trilock® STD Offset size 8
101204090	Trilock® STD Offset size 9
101204100	Trilock® STD Offset size 10
101204110	Trilock® STD Offset size 11
101204120	Trilock® STD Offset size 12
101214005	Trilock® LAT. size 0
101214010	Trilock® LAT. size 1
101214020	Trilock® LAT. size 2
101214030	Trilock® LAT. size 3
101214040	Trilock® LAT. size 4
101214050	Trilock® LAT. size 5
101214060	Trilock® LAT. size 6
101214070	Trilock® LAT. size 7
101214080	Trilock® LAT. size 8
101214090	Trilock® LAT. size 9
101214100	Trilock® LAT. size 10
101214110	Trilock® LAT. size 11
101214120	Trilock® LAT. size 12
121730038	Pinnacle Sector Gription 38mm
121730040	Pinnacle Sector Gription 40mm
121730042	Pinnacle Sector Gription 42mm
121730044	Pinnacle Sector Gription 44mm
121730046	Pinnacle Sector Gription 46mm
121732048	Pinnacle Sector Gription 48mm
121732050	Pinnacle Sector Gription 50mm

121732052	Pinnacle Sector Gription 52mm
121732054	Pinnacle Sector Gription 54mm
121732056	Pinnacle Sector Gription 56mm
121732058	Pinnacle Sector Gription 58mm
121732060	Pinnacle Sector Gription 60mm
121732062	Pinnacle Sector Gription 62mm
121732064	Pinnacle Sector Gription 64mm
121732066	Pinnacle Sector Gription 66mm
122128044	ALTRX Neutral 28x44mm
122128046	ALTRX Neutral 28x46mm
122128048	ALTRX Neutral 28x48mm
122128050	ALTRX Neutral 28x50mm
122128052	ALTRX Neutral 28x52mm
122128054	ALTRX Neutral 28x54mm
122132048	ALTRX Neutral 32x48mm
122132050	ALTRX Neutral 32x50mm
122132052	ALTRX Neutral 32x52mm
122132054	ALTRX Neutral 32x54mm
122132056	ALTRX Neutral 32x56mm
122132058	ALTRX Neutral 32x58mm
122136052	ALTRX Neutral 36x52mm
122136054	ALTRX Neutral 36x54mm
122136056	ALTRX Neutral 36x56mm
122136058	ALTRX Neutral 36x58mm
122136060	ALTRX Neutral 36x60mm
122136062	ALTRX Neutral 36x62mm
122136064	ALTRX Neutral 36x64mm
122136066	ALTRX Neutral 36x66mm
L92412	Stem diam. 10 lenght 225mm
L92414	Stem diam. 10 lenght 275mm
L92416	Stem diam. 10 lenght 325mm
L92418	Stem diam. 10 lenght 375mm
L92422	Stem diam. 12 lenght 225mm
L92424	Stem diam. 12 lenght 275mm
L92426	Stem diam. 12 lenght 325mm
L92428	Stem diam. 12 lenght 375mm
L92432	Stem diam. 14 lenght 225mm
L92434	Stem diam. 14 lenght 275mm
L92436	Stem diam. 14 lenght 325mm
L92438	Stem diam. 14 lenght 375mm
L92442	Stem diam. 16 lenght 225mm
L92444	Stem diam. 16 lenght 275mm
L92446	Stem diam. 16 lenght 325mm

L92448	Stem diam. 16 lenght 375mm
L92452	Stem diam. 18 lenght 225mm
L92454	Stem diam. 18 lenght 275mm
L92456	Stem diam. 18 lenght 325mm
L92458	Stem diam. 18 lenght 375mm
L92462	Stem diam. 20 lenght 225mm
L92464	Stem diam. 20 lenght 275mm
L92466	Stem diam. 20 lenght 325mm
L92468	Stem diam. 20 lenght 375mm
L92406	REEF™ modular neck without collar 25mm
L92409	REEF™ modular neck without collar 35mm
L92370	REEF™ Screw diam. 5mm lenght 20mm
L92372	REEF™ Screw diam. 5mm lenght 25mm
L92374	REEF™ Screw diam. 5mm lenght 30mm
L92376	REEF™ Screw diam. 5mm lenght 35mm
L92378	REEF™ Screw diam. 5mm lenght 40mm
L92380	REEF™ Screw diam. 5mm lenght 45mm
L92382	REEF™ Screw diam. 5mm lenght 50mm
L92384	REEF™ Screw diam. 5mm lenght 55mm
L92386	REEF™ Screw diam. 5mm lenght 60mm
L92388	REEF™ Screw diam. 5mm lenght 65mm
L92390	REEF™ Screw diam. 5mm lenght 70mm
L92392	REEF™ Screw diam. 5mm lenght 75mm
L92394	REEF™ Screw diam. 5mm lenght 80mm
L92412	Stem diam. 10 lenght 225mm
L92414	Stem diam. 10 lenght 275mm
L92416	Stem diam. 10 lenght 325mm
L92418	Stem diam. 10 lenght 375mm
L92422	Stem diam. 12 lenght 225mm
L92424	Stem diam. 12 lenght 275mm
L92426	Stem diam. 12 lenght 325mm
L92428	Stem diam. 12 lenght 375mm
L92432	Stem diam. 14 lenght 225mm
L92434	Stem diam. 14 lenght 275mm
L92436	Stem diam. 14 lenght 325mm
L92438	Stem diam. 14 lenght 375mm
L92442	Stem diam. 16 lenght 225mm
L92444	Stem diam. 16 lenght 275mm
L92446	Stem diam. 16 lenght 325mm
L92448	Stem diam. 16 lenght 375mm
L92452	Stem diam. 18 lenght 225mm
L92454	Stem diam. 18 lenght 275mm
L92456	Stem diam. 18 lenght 325mm

L92458	Stem diam. 18 length 375mm
L92462	Stem diam. 20 length 225mm
L92464	Stem diam. 20 length 275mm
L92466	Stem diam. 20 length 325mm
L92468	Stem diam. 20 length 375mm
L92406	REEF™ modular neck without collar 25mm
L92409	REEF™ modular neck without collar 35mm
L92370	REEF™ Screw diam. 5mm length 20mm
L92372	REEF™ Screw diam. 5mm length 25mm
L92374	REEF™ Screw diam. 5mm length 30mm
L92376	REEF™ Screw diam. 5mm length 35mm
L92378	REEF™ Screw diam. 5mm length 40mm
L92380	REEF™ Screw diam. 5mm length 45mm
L92382	REEF™ Screw diam. 5mm length 50mm
L92384	REEF™ Screw diam. 5mm length 55mm
L92386	REEF™ Screw diam. 5mm length 60mm
L92388	REEF™ Screw diam. 5mm length 65mm
L92390	REEF™ Screw diam. 5mm length 70mm
L92392	REEF™ Screw diam. 5mm length 75mm
L92394	REEF™ Screw diam. 5mm length 80mm
960007	P.F.C.® SIGMA CR non porous cemented femoral component left, size 1,5
960002	P.F.C.® SIGMA CR non porous cemented femoral component left, size 2
960008	P.F.C.® SIGMA CR non porous cemented femoral component left, size 2,5
960003	P.F.C.® SIGMA CR non porous cemented femoral component left, size 3
960004	P.F.C.® SIGMA CR non porous cemented femoral component left, size 4
196000400	P.F.C.® SIGMA CR non porous cemented femoral component left, size 4N
960005	P.F.C.® SIGMA CR non porous cemented femoral component left, size 5
960006	P.F.C.® SIGMA CR non porous cemented femoral component right, size 6
960017	P.F.C.® SIGMA CR non porous cemented femoral component right, size 1,5
960012	P.F.C.® SIGMA CR non porous cemented femoral component right, size 2
960018	P.F.C.® SIGMA CR non porous cemented femoral component right, size 2,5
960013	P.F.C.® SIGMA CR non porous cemented femoral component right, size 3
960014	P.F.C.® SIGMA CR non porous cemented femoral component right, size 4N
196001400	P.F.C.® SIGMA CR non porous cemented femoral component right, size 4N
960015	P.F.C.® SIGMA CR non porous cemented femoral component left, size 5
960016	P.F.C.® SIGMA CR non porous cemented femoral component right, size 6
940017	P.F.C.® SIGMA CR porous cemented femoral component left, size 1,5

940012	P.F.C.® SIGMA CR porous cemented femoral component left, size 2
940018	P.F.C.® SIGMA CR porous cemented femoral component left, size 2,5
940013	P.F.C.® SIGMA CR porous cemented femoral component left, size 3
940014	P.F.C.® SIGMA CR porous cemented femoral component left, size 4
940019	P.F.C.® SIGMA CR porous cemented femoral component left, size 4N
940015	P.F.C.® SIGMA CR porous cemented femoral component left, size 5
940016	P.F.C.® SIGMA CR porous cemented femoral component right, size 6
940027	P.F.C.® SIGMA CR porous cemented femoral component right, size 1,5
940022	P.F.C.® SIGMA CR porous cemented femoral component right, size 2
940028	P.F.C.® SIGMA CR porous cemented femoral component right, size 2,5
940023	P.F.C.® SIGMA CR porous cemented femoral component right, size 3
940024	P.F.C.® SIGMA CR porous cemented femoral component right, size 4N
940029	P.F.C.® SIGMA CR porous cemented femoral component right, size 4N
940025	P.F.C.® SIGMA CR porous cemented femoral component left, size 5
940026	P.F.C.® SIGMA CR porous cemented femoral component right, size 6
196040150	P.F.C.® SIGMA PS non porous cemented femoral component, left, size 1,5
196040200	P.F.C.® SIGMA PS non porous cemented femoral component, left, size 2
196040250	P.F.C.® SIGMA PS non porous cemented femoral component, left, size 2,5
196040300	P.F.C.® SIGMA PS non porous cemented femoral component, left, size 3
196040400	P.F.C.® SIGMA PS non porous cemented femoral component, left, size 4
196040450	P.F.C.® SIGMA PS non porous cemented femoral component, left, size 4N
196040500	P.F.C.® SIGMA PS non porous cemented femoral component, left, size 5
196040600	P.F.C.® SIGMA PS non porous cemented femoral component, right, size 6
196050150	P.F.C.® SIGMA PS non porous cemented femoral component, right, size 1,5
196050200	P.F.C.® SIGMA PS non porous cemented femoral component, right, size 2
196050250	P.F.C.® SIGMA PS non porous cemented femoral component, right, size 2,5
196050300	P.F.C.® SIGMA PS non porous cemented femoral component, right, size 3
196050400	P.F.C.® SIGMA PS non porous cemented femoral component, right, size 4
196050450	P.F.C.® SIGMA PS non porous cemented femoral component, right, size 4N
196050500	P.F.C.® SIGMA PS non porous cemented femoral component, left, size 5
196050600	P.F.C.® SIGMA PS non porous cemented femoral component, right, size 6
196020150	P.F.C.® SIGMA CR150 non porous, cemented, left, size 1,5
196020200	P.F.C.® SIGMA CR150 non porous, cemented, left, size 2

196020250	P.F.C.® SIGMA CR150 non porous, cemented, left, size 2,5
196020300	P.F.C.® SIGMA CR150 non porous, cemented, left, size 3
196020450	P.F.C.® SIGMA CR150 non porous, cemented, left, size 4N
196020400	P.F.C.® SIGMA CR150 non porous, cemented, left, size 4
196020500	P.F.C.® SIGMA CR150 non porous, cemented, left, size 5
196020600	P.F.C.® SIGMA CR150 non porous, cemented, left, size 6
196030150	P.F.C.® SIGMA CR150 non porous, cemented, right, size 1,5
196030200	P.F.C.® SIGMA CR150 non porous, cemented, right, size 2
196030250	P.F.C.® SIGMA CR150 non porous, cemented, right, size 2,5
196030300	P.F.C.® SIGMA CR150 non porous, cemented, right, size 3
196030450	P.F.C.® SIGMA CR150 non porous, cemented, right, size 4N
196030400	P.F.C.® SIGMA CR150 non porous, cemented, right, size 4
196030500	P.F.C.® SIGMA CR150 non porous, cemented, right, size 5
196030600	P.F.C.® SIGMA CR150 non porous, cemented, right, size 6
866027	P.F.C.® SIGMA Modular Ti Tibia tray, non porous, cemented, size 1,5
866022	P.F.C.® SIGMA Modular Ti Tibia tray, non porous, cemented, size 2
866029	P.F.C.® SIGMA Modular Ti Tibia tray, non porous, cemented, size 2,5
866023	P.F.C.® SIGMA Modular Ti Tibia tray, non porous, cemented, size 3
866024	P.F.C.® SIGMA Modular Ti Tibia tray, non porous, cemented, size 4
866025	P.F.C.® SIGMA Modular Ti Tibia tray, non porous, cemented, size 5
866026	P.F.C.® SIGMA Modular Ti Tibia tray, non porous, cemented, size 6
129433115	P.F.C.® SIGMA 'MBT Tibial tray, non porous, cemented size 1,5
129433120	P.F.C.® SIGMA 'MBT Tibial tray, non porous, cemented size 2
129433125	P.F.C.® SIGMA 'MBT Tibial tray, non porous, cemented size 2,5
129433130	P.F.C.® SIGMA 'MBT Tibial tray, non porous, cemented size 3
129433140	P.F.C.® SIGMA 'MBT Tibial tray, non porous, cemented size 4
129433150	P.F.C.® SIGMA 'MBT Tibial tray, non porous, cemented size 5
129433160	P.F.C.® SIGMA 'MBT Tibial tray, non porous, cemented size 6
860126	P.F.C.® SIGMA Modular plus Ti Tibia tray, non porous, cemented, size 1,5
860122	P.F.C.® SIGMA Modular plus Ti Tibia tray, non porous, cemented, size 2
860127	P.F.C.® SIGMA Modular plus Ti Tibia tray, non porous, cemented, size 2,5
860123	P.F.C.® SIGMA Modular plus Ti Tibia tray, non porous, cemented, size 3
860124	P.F.C.® SIGMA Modular plus Ti Tibia tray, non porous, cemented, size 4

860125	P.F.C.® SIGMA Modular plus Ti Tibia tray, non porous, cemented, size 5
158115000	P.F.C.® SIGMA CoCr Tibia tray non porous cemented, size 1,5
158120000	P.F.C.® SIGMA CoCr Tibia tray non porous cemented, size 2
158125000	P.F.C.® SIGMA CoCr Tibia tray non porous cemented, size 2,5
158130000	P.F.C.® SIGMA CoCr Tibia tray non porous cemented, size 3
158140000	P.F.C.® SIGMA CoCr Tibia tray non porous cemented, size 4
158150000	P.F.C.® SIGMA CoCr Tibia tray non porous cemented, size 5
158160000	P.F.C.® SIGMA CoCr Tibia tray non porous cemented, size 6
158110008	P.F.C.® SIGMA CR CVD GVF tibial insert, size 1,5 8mm
158110010	P.F.C.® SIGMA CR CVD GVF tibial insert, size 1,5 10mm
158110012	P.F.C.® SIGMA CR CVD GVF tibial insert, size 1,5 12,5mm
158110015	P.F.C.® SIGMA CR CVD GVF tibial insert, size 1,5 15mm
158110017	P.F.C.® SIGMA CR CVD GVF tibial insert, size 1,5 17,5mm
158110020	P.F.C.® SIGMA CR CVD GVF tibial insert, size 1,5 20mm
158111008	P.F.C.® SIGMA CR CVD GVF tibial insert, size 2 8mm
158111010	P.F.C.® SIGMA CR CVD GVF tibial insert, size 2 10mm
158111012	P.F.C.® SIGMA CR CVD GVF tibial insert, size 2 12,5mm
158111015	P.F.C.® SIGMA CR CVD GVF tibial insert, size 2 15mm
158111017	P.F.C.® SIGMA CR CVD GVF tibial insert, size 2 17,5mm
158111020	P.F.C.® SIGMA CR CVD GVF tibial insert, size 2 20mm
158112008	P.F.C.® SIGMA CR CVD GVF tibial insert, size 2,5 8mm
158112010	P.F.C.® SIGMA CR CVD GVF tibial insert, size 2,5 10mm
158112012	P.F.C.® SIGMA CR CVD GVF tibial insert, size 2,5 12,5mm
158112015	P.F.C.® SIGMA CR CVD GVF tibial insert, size 2,5 15mm
158112017	P.F.C.® SIGMA CR CVD GVF tibial insert, size 2,5 17,5mm
158112020	P.F.C.® SIGMA CR CVD GVF tibial insert, size 2,5 20mm
158113008	P.F.C.® SIGMA CR CVD GVF tibial insert, size 3 8mm
158113010	P.F.C.® SIGMA CR CVD GVF tibial insert, size 3 10mm
158113012	P.F.C.® SIGMA CR CVD GVF tibial insert, size 3 12,5mm
158113015	P.F.C.® SIGMA CR CVD GVF tibial insert, size 3 15mm
158113017	P.F.C.® SIGMA CR CVD GVF tibial insert, size 3 17,5mm
158113020	P.F.C.® SIGMA CR CVD GVF tibial insert, size 3 20mm
158114008	P.F.C.® SIGMA CR CVD GVF tibial insert, size 4 8mm
158114010	P.F.C.® SIGMA CR CVD GVF tibial insert, size 4 10mm

158114012	P.F.C.® SIGMA CR CVD GVF tibial insert, size 4 12,5mm
158114015	P.F.C.® SIGMA CR CVD GVF tibial insert, size 4 15mm
158114017	P.F.C.® SIGMA CR CVD GVF tibial insert, size 4 17,5mm
158114020	P.F.C.® SIGMA CR CVD GVF tibial insert, size 4 20mm
158115008	P.F.C.® SIGMA CR CVD GVF tibial insert, size 5 8mm
158115010	P.F.C.® SIGMA CR CVD GVF tibial insert, size 5 10mm
158115012	P.F.C.® SIGMA CR CVD GVF tibial insert, size 5 12,5mm
158115015	P.F.C.® SIGMA CR CVD GVF tibial insert, size 5 15mm
158115017	P.F.C.® SIGMA CR CVD GVF tibial insert, size 5 17,5mm
158115020	P.F.C.® SIGMA CR CVD GVF tibial insert, size 5 20mm
158116008	P.F.C.® SIGMA CR CVD GVF tibial insert, size 6 8mm
158116010	P.F.C.® SIGMA CR CVD GVF tibial insert, size 6 10mm
158116012	P.F.C.® SIGMA CR CVD GVF tibial insert, size 6 12,5mm
158116015	P.F.C.® SIGMA CR CVD GVF tibial insert, size 6 15mm
158116017	P.F.C.® SIGMA CR CVD GVF tibial insert, size 6 17,5mm
158116020	P.F.C.® SIGMA CR CVD GVF tibial insert, size 6 20mm
972001	P.F.C.® SIGMA CR CVD RP GVF tibial insert, size 1,5 10mm
972002	P.F.C.® SIGMA CR CVD RP GVF tibial insert, size 1,5 12,5mm
972003	P.F.C.® SIGMA CR CVD RP GVF tibial insert, size 1,5 15mm
972004	P.F.C.® SIGMA CR CVD RP GVF tibial insert, size 1,5 17,5mm
962011	P.F.C.® SIGMA CR CVD RP GVF tibial insert, size 2 10mm
962012	P.F.C.® SIGMA CR CVD RP GVF tibial insert, size 2 12,5mm
962013	P.F.C.® SIGMA CR CVD RP GVF tibial insert, size 2 15mm
962014	P.F.C.® SIGMA CR CVD RP GVF tibial insert, size 2 17,5mm
962021	P.F.C.® SIGMA CR CVD RP GVF tibial insert, size 2,5 10mm
962022	P.F.C.® SIGMA CR CVD RP GVF tibial insert, size 2,5 12,5mm
962023	P.F.C.® SIGMA CR CVD RP GVF tibial insert, size 2,5 15mm
962024	P.F.C.® SIGMA CR CVD RP GVF tibial insert, size 2,5 17mm
962031	P.F.C.® SIGMA CR CVD RP GVF tibial insert, size 3 10mm
962032	P.F.C.® SIGMA CR CVD RP GVF tibial insert, size 3 12,5mm
962033	P.F.C.® SIGMA CR CVD RP GVF tibial insert, size 3 15mm
962034	P.F.C.® SIGMA CR CVD RP GVF tibial insert, size 3 17,5mm
962041	P.F.C.® SIGMA CR CVD RP GVF tibial insert, size 4 10mm

962042	P.F.C.® SIGMA CR CVD RP GVF tibial insert, size 4 12,5mm
962043	P.F.C.® SIGMA CR CVD RP GVF tibial insert, size 4 15mm
962044	P.F.C.® SIGMA CR CVD RP GVF tibial insert, size 4 17,5mm
962051	P.F.C.® SIGMA CR CVD RP GVF tibial insert, size 5 10mm
962052	P.F.C.® SIGMA CR CVD RP GVF tibial insert, size 5 12,5mm
962053	P.F.C.® SIGMA CR CVD RP GVF tibial insert, size 5 15mm
962054	P.F.C.® SIGMA CR CVD RP GVF tibial insert, size 5 17mm
962061	P.F.C.® SIGMA CR CVD RP GVF tibial insert, size 6 10mm
962062	P.F.C.® SIGMA CR CVD RP GVF tibial insert, size 6 12,5mm
962063	P.F.C.® SIGMA CR CVD RP GVF tibial insert, size 6 15mm
962064	P.F.C.® SIGMA CR CVD RP GVF tibial insert, size 6 17mm
158120008	P.F.C.® SIGMA STB GVF tibial insert, size 1.5 8mm
158120010	P.F.C.® SIGMA STB GVF tibial insert, size 1.5 10mm
158120012	P.F.C.® SIGMA STB GVF tibial insert, size 1.5 12,5mm
158120015	P.F.C.® SIGMA STB GVF tibial insert, size 1.5 15mm
158120017	P.F.C.® SIGMA STB GVF tibial insert, size 1.5 17,5mm
158121008	P.F.C.® SIGMA STB GVF tibial insert, size 2 8mm
158121010	P.F.C.® SIGMA STB GVF tibial insert, size 2 10mm
158121012	P.F.C.® SIGMA STB GVF tibial insert, size 2 12,5mm
158121015	P.F.C.® SIGMA STB GVF tibial insert, size 2 15mm
158121017	P.F.C.® SIGMA STB GVF tibial insert, size 2 17,5mm
158121020	P.F.C.® SIGMA STB GVF tibial insert, size 2 20mm
158121022	P.F.C.® SIGMA STB GVF tibial insert, size 2 22,5mm
158121025	P.F.C.® SIGMA STB GVF tibial insert, size 2 25mm
158122008	P.F.C.® SIGMA STB GVF tibial insert, size 2,5 8mm
158122010	P.F.C.® SIGMA STB GVF tibial insert, size 2,5 10mm
158122012	P.F.C.® SIGMA STB GVF tibial insert, size 2,5 12,5mm
158122015	P.F.C.® SIGMA STB GVF tibial insert, size 2,5 15mm
158122017	P.F.C.® SIGMA STB GVF tibial insert, size 2,5 17,5mm
158122020	P.F.C.® SIGMA STB GVF tibial insert, size 2,5 20mm
158122022	P.F.C.® SIGMA STB GVF tibial insert, size 2,5 22,5mm
158122025	P.F.C.® SIGMA STB GVF tibial insert, size 2,5 25mm
158123008	P.F.C.® SIGMA STB GVF tibial insert, size 3 8mm
158123010	P.F.C.® SIGMA STB GVF tibial insert, size 3 10mm
158123012	P.F.C.® SIGMA STB GVF tibial insert, size 3 12,5mm
158123015	P.F.C.® SIGMA STB GVF tibial insert, size 3 15mm
158123017	P.F.C.® SIGMA STB GVF tibial insert, size 3 17,5mm
158123020	P.F.C.® SIGMA STB GVF tibial insert, size 3 20mm

158123022	P.F.C.® SIGMA STB GVF tibial insert, size 3 22,5mm
158123025	P.F.C.® SIGMA STB GVF tibial insert, size 3 25mm
158124008	P.F.C.® SIGMA STB GVF tibial insert, size 4 8mm
158124010	P.F.C.® SIGMA STB GVF tibial insert, size 4 10mm
158124012	P.F.C.® SIGMA STB GVF tibial insert, size 4 12,5mm
158124015	P.F.C.® SIGMA STB GVF tibial insert, size 4 15mm
158124017	P.F.C.® SIGMA STB GVF tibial insert, size 4 17,5mm
158124020	P.F.C.® SIGMA STB GVF tibial insert, size 4 20mm
158124022	P.F.C.® SIGMA STB GVF tibial insert, size 4 22,5mm
158124025	P.F.C.® SIGMA STB GVF tibial insert, size 4 25mm
158125008	P.F.C.® SIGMA STB GVF tibial insert, size 5 8mm
158125010	P.F.C.® SIGMA STB GVF tibial insert, size 5 10mm
158125012	P.F.C.® SIGMA STB GVF tibial insert, size 5 12,5mm
158125015	P.F.C.® SIGMA STB GVF tibial insert, size 5 15mm
158125017	P.F.C.® SIGMA STB GVF tibial insert, size 5 17,5mm
158125020	P.F.C.® SIGMA STB GVF tibial insert, size 5 20mm
158125022	P.F.C.® SIGMA STB GVF tibial insert, size 5 22,5mm
158125025	P.F.C.® SIGMA STB GVF tibial insert, size 5 25mm
158126008	P.F.C.® SIGMA STB GVF tibial insert, size 6 8mm
158126010	P.F.C.® SIGMA STB GVF tibial insert, size 6 10mm
158126012	P.F.C.® SIGMA STB GVF tibial insert, size 6 12,5mm
158126015	P.F.C.® SIGMA STB GVF tibial insert, size 6 15mm
158126017	P.F.C.® SIGMA STB GVF tibial insert, size 6 17,5mm
158126020	P.F.C.® SIGMA STB GVF tibial insert, size 6 20mm
158126022	P.F.C.® SIGMA STB GVF tibial insert, size 6 22,5mm
158126025	P.F.C.® SIGMA STB GVF tibial insert, size 6 25mm
972101	P.F.C.® SIGMA PS STB RP GVF tibial insert, size 1,5 10mm
972102	P.F.C.® SIGMA PS STB RP GVF tibial insert, size 1,5 12,5mm
972103	P.F.C.® SIGMA PS STB RP GVF tibial insert, size 1,5 15mm
972104	P.F.C.® SIGMA PS STB RP GVF tibial insert, size 1,5 17,5mm
972105	P.F.C.® SIGMA PS STB RP GVF tibial insert, size 1,5 20mm
962111	P.F.C.® SIGMA PS STB RP GVF tibial insert, size 2 10mm
962112	P.F.C.® SIGMA PS STB RP GVF tibial insert, size 2 12,5mm
962113	P.F.C.® SIGMA PS STB RP GVF tibial insert, size 2 15mm
962114	P.F.C.® SIGMA PS STB RP GVF tibial insert, size 2 17,5mm
962115	P.F.C.® SIGMA PS STB RP GVF tibial insert, size 2 20mm
962116	P.F.C.® SIGMA PS STB RP GVF tibial insert, size 2 22,5mm
962117	P.F.C.® SIGMA PS STB RP GVF tibial insert, size 2 25mm
962121	P.F.C.® SIGMA PS STB RP GVF tibial insert, size 2,5 10mm

962122	P.F.C.® SIGMA PS STB RP GVF tibial insert, size 2,5 12,5mm
962123	P.F.C.® SIGMA PS STB RP GVF tibial insert, size 2,5 15mm
962124	P.F.C.® SIGMA PS STB RP GVF tibial insert, size 2,5 17mm
962125	P.F.C.® SIGMA PS STB RP GVF tibial insert, size 2,5 20mm
962126	P.F.C.® SIGMA PS STB RP GVF tibial insert, size 2,5 22,5mm
962127	P.F.C.® SIGMA PS STB RP GVF tibial insert, size 2,5 25mm
962131	P.F.C.® SIGMA PS STB RP GVF tibial insert, size 3 10mm
962132	P.F.C.® SIGMA PS STB RP GVF tibial insert, size 3 12,5mm
962133	P.F.C.® SIGMA PS STB RP GVF tibial insert, size 3 15mm
962134	P.F.C.® SIGMA PS STB RP GVF tibial insert, size 3 17,5mm
962135	P.F.C.® SIGMA PS STB RP GVF tibial insert, size 3 20mm
962136	P.F.C.® SIGMA PS STB RP GVF tibial insert, size 3 22,5mm
962137	P.F.C.® SIGMA PS STB RP GVF tibial insert, size 3 25mm
962141	P.F.C.® SIGMA PS STB RP GVF tibial insert, size 4 10mm
962142	P.F.C.® SIGMA PS STB RP GVF tibial insert, size 4 12,5mm
962143	P.F.C.® SIGMA PS STB RP GVF tibial insert, size 4 15mm
962144	P.F.C.® SIGMA PS STB RP GVF tibial insert, size 4 17,5mm
962145	P.F.C.® SIGMA PS STB RP GVF tibial insert, size 4 20mm
962146	P.F.C.® SIGMA PS STB RP GVF tibial insert, size 4 22,5mm
962147	P.F.C.® SIGMA PS STB RP GVF tibial insert, size 4 25mm
962151	P.F.C.® SIGMA PS STB RP GVF tibial insert, size 5 10mm
962152	P.F.C.® SIGMA PS STB RP GVF tibial insert, size 5 12,5mm
962153	P.F.C.® SIGMA PS STB RP GVF tibial insert, size 5 15mm
962154	P.F.C.® SIGMA PS STB RP GVF tibial insert, size 5 17,5mm
962155	P.F.C.® SIGMA PS STB RP GVF tibial insert, size 5 20mm
962156	P.F.C.® SIGMA PS STB RP GVF tibial insert, size 5 22,5mm
962157	P.F.C.® SIGMA PS STB RP GVF tibial insert, size 5 25mm
962161	P.F.C.® SIGMA PS STB RP GVF tibial insert, size 6 10mm
962162	P.F.C.® SIGMA PS STB RP GVF tibial insert, size 6 12,5mm
962163	P.F.C.® SIGMA PS STB RP GVF tibial insert, size 6 15mm
962164	P.F.C.® SIGMA PS STB RP GVF tibial insert, size 6 17mm
962165	P.F.C.® SIGMA PS STB RP GVF tibial insert, size 6 20mm
962166	P.F.C.® SIGMA PS STB RP GVF tibial insert, size 6 22,5mm

962167	P.F.C.® SIGMA PS STB RP GVF tibial insert, size 6 25mm
960100	P.F.C.® SIGMA Oval Dome patellar prosthesis 3 peg 32mm
960101	P.F.C.® SIGMA Oval Dome patellar prosthesis 3 peg 35mm
960102	P.F.C.® SIGMA Oval Dome patellar prosthesis 3 peg 38mm
960103	P.F.C.® SIGMA Oval Dome patellar prosthesis 3 peg 41mm
960085	P.F.C.® SIGMA TC3 femur, left, size 1,5
960081	P.F.C.® SIGMA TC3 femur, left, size 2
960086	P.F.C.® SIGMA TC3 femur, left, size 2,5
960082	P.F.C.® SIGMA TC3 femur, left, size 3
960083	P.F.C.® SIGMA TC3 femur, left, size 4
960084	P.F.C.® SIGMA TC3 femur, left, size 5
960091	P.F.C.® SIGMA TC3 femur, right, size 1,5
960087	P.F.C.® SIGMA TC3 femur, right, size 2
960092	P.F.C.® SIGMA TC3 femur, right, size 2,5
960088	P.F.C.® SIGMA TC3 femur, right, size 3
960089	P.F.C.® SIGMA TC3 femur, right, size 4
960090	P.F.C.® SIGMA TC3 femur, right, size 5
867410	P.F.C.® Sigma fluted stem 10x75 mm
867412	P.F.C.® Sigma fluted stem 12x75 mm
867414	P.F.C.® Sigma fluted stem 14x75 mm
867416	P.F.C.® Sigma fluted stem 16x75 mm
867418	P.F.C.® Sigma fluted stem 18x75 mm
867419	P.F.C.® Sigma fluted stem 20x75 mm
867420	P.F.C.® Sigma fluted stem 22x75 mm
867421	P.F.C.® Sigma fluted stem 24x75 mm
867424	P.F.C.® Sigma fluted stem 10x115 mm
867426	P.F.C.® Sigma fluted stem 12x115 mm
867428	P.F.C.® Sigma fluted stem 14x115 mm
867430	P.F.C.® Sigma fluted stem 16x115 mm
867432	P.F.C.® Sigma fluted stem 18x115 mm
867433	P.F.C.® Sigma fluted stem 20x115 mm
867434	P.F.C.® Sigma fluted stem 22x115 mm
867435	P.F.C.® Sigma fluted stem 24x115 mm
867438	P.F.C.® Sigma fluted stem 10x150 mm
867440	P.F.C.® Sigma fluted stem 12x150 mm
867442	P.F.C.® Sigma fluted stem 14x150 mm
867444	P.F.C.® Sigma fluted stem 16x150 mm
867446	P.F.C.® Sigma fluted stem 18x150 mm
867447	P.F.C.® Sigma fluted stem 20x150 mm
867448	P.F.C.® Sigma fluted stem 22x150 mm
867449	P.F.C.® Sigma fluted stem 24x150 mm

866401	P.F.C.® SIGMA cemented stem 13x30 mm
866402	P.F.C.® SIGMA cemented stem 13x60 mm
866403	P.F.C.® SIGMA cemented stem 15x30 mm
866404	P.F.C.® SIGMA cemented stem 15x60 mm
866468	P.F.C.® SIGMA cemented stem 13x90 mm
866469	P.F.C.® SIGMA cemented stem 15x90 mm
129435115	P.F.C.® Sigma MBT revision cemented tibia, size 1,5
129435120	P.F.C.® Sigma MBT revision cemented tibia, size 2
129435125	P.F.C.® Sigma MBT revision cemented tibia, size 2,5
129435130	P.F.C.® Sigma MBT revision cemented tibia, size 3
129435140	P.F.C.® Sigma MBT revision cemented tibia, size 4
129435150	P.F.C.® Sigma MBT revision cemented tibia, size 5
129435160	P.F.C.® Sigma MBT revision cemented tibia, size 6
960510	P.F.C.® Sigma TC3 insert 10mm, velikost 2
960511	P.F.C.® Sigma TC3 insert 12mm, velikost 2
960512	P.F.C.® Sigma TC3 insert 15mm, velikost 2
960513	P.F.C.® Sigma TC3 insert 17mm, velikost 2
960514	P.F.C.® Sigma TC3 insert 20mm, velikost 2
960515	P.F.C.® Sigma TC3 insert 22mm, velikost 2
960516	P.F.C.® Sigma TC3 insert 25mm, velikost 2
960517	P.F.C.® Sigma TC3 insert 30mm, velikost 2
960530	P.F.C.® Sigma TC3 insert 10mm, velikost 3
960531	P.F.C.® Sigma TC3 insert 12mm, velikost 3
960532	P.F.C.® Sigma TC3 insert 15mm, velikost 3
960533	P.F.C.® Sigma TC3 insert 17mm, velikost 3
960534	P.F.C.® Sigma TC3 insert 20mm, velikost 3
960535	P.F.C.® Sigma TC3 insert 22mm, velikost 3
960536	P.F.C.® Sigma TC3 insert 25mm, velikost 3
960537	P.F.C.® Sigma TC3 insert 30mm, velikost 3
960540	P.F.C.® Sigma TC3 insert 10mm, velikost 4
960541	P.F.C.® Sigma TC3 insert 12mm, velikost 4
960542	P.F.C.® Sigma TC3 insert 15mm, velikost 4
960543	P.F.C.® Sigma TC3 insert 17mm, velikost 4
960544	P.F.C.® Sigma TC3 insert 20mm, velikost 4
960545	P.F.C.® Sigma TC3 insert 22mm, velikost 4
960546	P.F.C.® Sigma TC3 insert 25mm, velikost 4
960550	P.F.C.® Sigma TC3 insert 10mm, velikost 5
960551	P.F.C.® Sigma TC3 insert 12mm, velikost 5
960552	P.F.C.® Sigma TC3 insert 15mm, velikost 5
960553	P.F.C.® Sigma TC3 insert 17mm, velikost 5
960554	P.F.C.® Sigma TC3 insert 20mm, velikost 5
960555	P.F.C.® Sigma TC3 insert 22mm, velikost 5

960556	P.F.C.® Sigma TC3 insert 25mm, velikost 5
960557	P.F.C.® Sigma TC3 insert 30mm, velikost 5
962311	P.F.C.® Sigma TC3 RP insert 1.5 -10mm
962312	P.F.C.® Sigma TC3 RP insert 1.5 -12.5mm
962313	P.F.C.® Sigma TC3 RP insert 1.5 -15mm
962314	P.F.C.® Sigma TC3 RP insert 1.5 -17.5mm
962315	P.F.C.® Sigma TC3 RP insert 1.5 -20mm
962316	P.F.C.® Sigma TC3 RP insert 1.5 -22.5mm
962317	P.F.C.® Sigma TC3 RP insert 1.5 -25mm
962318	P.F.C.® Sigma TC3 RP insert 1.5 -30mm
962321	P.F.C.® Sigma TC3 RP insert 2 -10mm
962322	P.F.C.® Sigma TC3 RP insert 2 -12.5mm
962323	P.F.C.® Sigma TC3 RP insert 2- 15mm
962324	P.F.C.® Sigma TC3 RP insert 2-17.5mm
962325	P.F.C.® Sigma TC3 RP insert 2-20mm
962326	P.F.C.® Sigma TC3 RP insert 2-22.5mm
962327	P.F.C.® Sigma TC3 RP insert 2-25mm
962328	P.F.C.® Sigma TC3 RP insert 2-30mm
962331	P.F.C.® Sigma TC3 RP insert 2.5 -10mm
962332	P.F.C.® Sigma TC3 RP insert 2.5- 12.5mm
962333	P.F.C.® Sigma TC3 RP insert 2.5- 15mm
962334	P.F.C.® Sigma TC3 RP insert 2.5- 17.5mm
962335	P.F.C.® Sigma TC3 RP insert 2.5- 20mm
962336	P.F.C.® Sigma TC3 RP insert 2.5-22.5mm
962337	P.F.C.® Sigma TC3 RP insert 2.5- 25mm
962338	P.F.C.® Sigma TC3 RP insert 2.5- 30mm
962341	P.F.C.® Sigma TC3 RP insert 3-10mm
962342	P.F.C.® Sigma TC3 RP insert 3- 12.5mm
962343	P.F.C.® Sigma TC3 RP insert 3- 15mm
962344	P.F.C.® Sigma TC3 RP insert 3-17.5mm
962345	P.F.C.® Sigma TC3 RP insert 3- 20mm
962346	P.F.C.® Sigma TC3 RP insert 3- 22.5mm
962347	P.F.C.® Sigma TC3 RP insert 3- 25.0mm
962348	P.F.C.® Sigma TC3 RP insert 3-30.0mm
962351	P.F.C.® Sigma TC3 RP insert 4- 10mm
962352	P.F.C.® Sigma TC3 RP insert 4- 12.5mm
962353	P.F.C.® Sigma TC3 RP insert 4- 15mm
962354	P.F.C.® Sigma TC3 RP insert 4-17.5mm
962355	P.F.C.® Sigma TC3 RP insert 4- 20mm
962356	P.F.C.® Sigma TC3 RP insert 4-22.5mm
962357	P.F.C.® Sigma TC3 RP insert 4-25mm
962358	P.F.C.® Sigma TC3 RP insert 4- 30mm

962361	P.F.C.® Sigma TC3 RP insert 5- 10mm
962362	P.F.C.® Sigma TC3 RP insert 5-12.5mm
962363	P.F.C.® Sigma TC3 RP insert 5- 15mm
962364	P.F.C.® Sigma TC3 RP insert 5-17.5mm
962365	P.F.C.® Sigma TC3 RP insert 5-20mm
962366	P.F.C.® Sigma TC3 RP insert 5-22.5mm
962367	P.F.C.® Sigma TC3 RP insert 5-25mm
962368	P.F.C.® Sigma TC3 RP insert 5-30mm
129456115	P.F.C.® Sigma MBT REVISION STEP WEDGE SZ1.5 5MM
129456116	P.F.C.® Sigma MBT REVISION STEP WEDGE SZ1.5 10MM
129456117	P.F.C.® Sigma MBT REVISION STEP WEDGE SZ1.5 15MM
129456120	P.F.C.® Sigma MBT REVISION STEP WEDGE SZ2 5MM
129456121	P.F.C.® Sigma MBT REVISION STEP WEDGE SZ2 10MM
129456122	P.F.C.® Sigma MBT REVISION STEP WEDGE SZ2 15MM
129456125	P.F.C.® Sigma MBT REVISION STEP WEDGE SZ2.5 5MM
129456126	P.F.C.® Sigma MBT REVISION STEP WEDGE SZ2.5 10MM
129456127	P.F.C.® Sigma MBT REVISION STEP WEDGE SZ2.5 15MM
129456130	P.F.C.® Sigma MBT REVISION STEP WEDGE SZ3 5MM
129456131	P.F.C.® Sigma MBT REVISION STEP WEDGE SZ3 10MM
129456132	P.F.C.® Sigma MBT REVISION STEP WEDGE SZ3 15MM
129456135	P.F.C.® Sigma MBT REVISION STEP WEDGE SZ4 5MM
129456136	P.F.C.® Sigma MBT REVISION STEP WEDGE SZ4 10MM
129456137	P.F.C.® Sigma MBT REVISION STEP WEDGE SZ4 15MM
129456140	P.F.C.® Sigma MBT REVISION STEP WEDGE SZ5 5MM
129456141	P.F.C.® Sigma MBT REVISION STEP WEDGE SZ5 10MM
129456142	P.F.C.® Sigma MBT REVISION STEP WEDGE SZ5 15MM
960820	P.F.C.® SIGMA distal femoral augment 4 mm, levi št. 2
960821	P.F.C.® SIGMA distal femoral augment 4 mm, desni št. 2
960822	P.F.C.® SIGMA distal femoral augment 8 mm, levi št. 2
960823	P.F.C.® SIGMA distal femoral augment 8 mm, desni št. 2
960826	P.F.C.® SIGMA posterior femoral augment 4 mm, št. 2
960828	P.F.C.® SIGMA posterior femoral augment 8 mm, št. 2
960829	P.F.C.® SIGMA distal femoral augment 12 mm, levi št. 2
960830	P.F.C.® SIGMA distal femoral augment 12 mm, desni št. 2
960831	P.F.C.® SIGMA distal femoral augment 16 mm, levi št. 2

960832	P.F.C.® SIGMA distal femoral augment 16 mm, desni št. 2
960860	P.F.C.® SIGMA distal femoral augment 4 mm, levi št. 3
960861	P.F.C.® SIGMA distal femoral augment 4 mm, desni št. 3
960862	P.F.C.® SIGMA distal femoral augment 8 mm, levi št. 3
960863	P.F.C.® SIGMA distal femoral augment 8 mm, desni št. 3
960866	P.F.C.® SIGMA posterior femoral augment 4 mm, št. 3
960868	P.F.C.® SIGMA posterior femoral augment 8 mm, št. 3
960869	P.F.C.® SIGMA distal femoral augment 12 mm, levi št. 3
960870	P.F.C.® SIGMA distal femoral augment 12 mm, desni št. 3
960871	P.F.C.® SIGMA distal femoral augment 16 mm, levi št. 3
960872	P.F.C.® SIGMA distal femoral augment 16 mm, desni št. 3
960880	P.F.C.® SIGMA distal femoral augment 4 mm, levi št. 4
960881	P.F.C.® SIGMA distal femoral augment 4 mm, desni št. 4
960882	P.F.C.® SIGMA distal femoral augment 8 mm, levi št. 4
960883	P.F.C.® SIGMA distal femoral augment 8 mm, desni št. 4
960886	P.F.C.® SIGMA posterior femoral augment 4 mm, št. 4
960888	P.F.C.® SIGMA posterior femoral augment 8 mm, št. 4
960889	P.F.C.® SIGMA distal femoral augment 12 mm, levi št. 4
960890	P.F.C.® SIGMA distal femoral augment 12 mm, desni št. 4
960891	P.F.C.® SIGMA distal femoral augment 16 mm, levi št. 4
960892	P.F.C.® SIGMA distal femoral augment 16 mm, desni št. 4
960900	P.F.C.® SIGMA distal femoral augment 4 mm, levi št. 5
960901	P.F.C.® SIGMA distal femoral augment 4 mm, desni št. 5
960902	P.F.C.® SIGMA distal femoral augment 8 mm, levi št. 5
960903	P.F.C.® SIGMA distal femoral augment 8 mm, desni št. 5
960906	P.F.C.® SIGMA posterior femoral augment 4 mm, št. 5
960908	P.F.C.® SIGMA posterior femoral augment 8 mm, št. 5
960909	P.F.C.® SIGMA distal femoral augment 12 mm, levi št. 5
960910	P.F.C.® SIGMA distal femoral augment 12 mm, desni št. 5
960911	P.F.C.® SIGMA distal femoral augment 16 mm, levi št. 5
960912	P.F.C.® SIGMA distal femoral augment 16 mm, desni št. 5
129453205	P.F.C.® SIGMA Femoral Sleeve 20mm
129453215	P.F.C.® SIGMA Femoral Sleeve 31mm
129453225	P.F.C.® SIGMA Femoral Sleeve 34mm

129453235	P.F.C.® SIGMA Femoral Sleeve 40mm
129453245	P.F.C.® SIGMA Femoral Sleeve 46mm
129454140	P.F.C.® SIGMA MBT Tibial Sleeve 29MM
129454100	P.F.C.® SIGMA MBT Tibial Sleeve 37MM
129454110	P.F.C.® SIGMA MBT Tibial Sleeve 45MM
129454120	P.F.C.® SIGMA MBT Tibial Sleeve 53MM
129454130	P.F.C.® SIGMA MBT Tibial Sleeve 61MM
960783	P.F.C.® SIGMA Neutral bolt for femoral adapter
960784	P.F.C.® SIGMA Offset bolt for femoral adapter
960781	P.F.C.® SIGMA Femoral adapter 5 degrees
960782	P.F.C.® SIGMA Femoral adapter 7 degrees
150400103	ATTUNE CR FEM LT SZ 3 CEM
150400104	ATTUNE CR FEM LT SZ 4 CEM
150400105	ATTUNE CR FEM LT SZ 5 CEM
150400106	ATTUNE CR FEM LT SZ 6 CEM
150400107	ATTUNE CR FEM LT SZ 7 CEM
150400108	ATTUNE CR FEM LT SZ 8 CEM
150400123	ATTUNE CR FEM LT SZ 3 NAR CEM
150400124	ATTUNE CR FEM LT SZ 4 NAR CEM
150400125	ATTUNE CR FEM LT SZ 5 NAR CEM
150400126	ATTUNE CR FEM LT SZ 6 NAR CEM
150400101	ATTUNE CR FEM LT SZ 1 CEM
150400102	ATTUNE CR FEM LT SZ 2 CEM
150400109	ATTUNE CR FEM LT SZ 9 CEM
150400110	ATTUNE CR FEM LT SZ 10 CEM
150400203	ATTUNE CR FEM RT SZ 3 CEM
150400204	ATTUNE CR FEM RT SZ 4 CEM
150400205	ATTUNE CR FEM RT SZ 5 CEM
150400206	ATTUNE CR FEM RT SZ 6 CEM
150400207	ATTUNE CR FEM RT SZ 7 CEM
150400208	ATTUNE CR FEM RT SZ 8 CEM
150400223	ATTUNE CR FEM RT SZ 3 NAR CEM
150400224	ATTUNE CR FEM RT SZ 4 NAR CEM
150400225	ATTUNE CR FEM RT SZ 5 NAR CEM
150400226	ATTUNE CR FEM RT SZ 6 NAR CEM
150400201	ATTUNE CR FEM RT SZ 1 CEM
150400202	ATTUNE CR FEM RT SZ 2 CEM
150400209	ATTUNE CR FEM RT SZ 9 CEM
150400210	ATTUNE CR FEM RT SZ 10 CEM
150410103	ATTUNE PS FEM LT SZ 3 CEM
150410104	ATTUNE PS FEM LT SZ 4 CEM
150410105	ATTUNE PS FEM LT SZ 5 CEM
150410106	ATTUNE PS FEM LT SZ 6 CEM

150410107	ATTUNE PS FEM LT SZ 7 CEM
150410108	ATTUNE PS FEM LT SZ 8 CEM
150410123	ATTUNE PS FEM LT SZ 3 NAR CEM
150410124	ATTUNE PS FEM LT SZ 4 NAR CEM
150410125	ATTUNE PS FEM LT SZ 5 NAR CEM
150410126	ATTUNE PS FEM LT SZ 6 NAR CEM
150410101	ATTUNE PS FEM LT SZ 1 CEM
150410102	ATTUNE PS FEM LT SZ 2 CEM
150410109	ATTUNE PS FEM LT SZ 9 CEM
150410110	ATTUNE PS FEM LT SZ 10 CEM
150410203	ATTUNE PS FEM RT SZ 3 CEM
150410204	ATTUNE PS FEM RT SZ 4 CEM
150410205	ATTUNE PS FEM RT SZ 5 CEM
150410206	ATTUNE PS FEM RT SZ 6 CEM
150410207	ATTUNE PS FEM RT SZ 7 CEM
150410208	ATTUNE PS FEM RT SZ 8 CEM
150410223	ATTUNE PS FEM RT SZ 3 NAR CEM
150410224	ATTUNE PS FEM RT SZ 4 NAR CEM
150410225	ATTUNE PS FEM RT SZ 5 NAR CEM
150410226	ATTUNE PS FEM RT SZ 6 NAR CEM
150410201	ATTUNE PS FEM RT SZ 9 CEM
150410202	ATTUNE PS FEM RT SZ 10 CEM
150401103	ATT CR FEM CLS SZ 3 LEFT
150401104	ATT CR FEM CLS SZ 4 LEFT
150401105	ATT CR FEM CLS SZ 5 LEFT
150401106	ATT CR FEM CLS SZ 6 LEFT
150401107	ATT CR FEM CLS SZ 7 LEFT
150401108	ATT CR FEM CLS SZ 8 LEFT
150401123	ATT CR FEM CLS SZ 3 N LEFT
150401124	ATT CR FEM CLS SZ 4 N LEFT
150401125	ATT CR FEM CLS SZ 5 N LEFT
150401126	ATT CR FEM CLS SZ 6 N LEFT
150401203	ATT CR FEM CLS SZ 3 RIGHT
150401204	ATT CR FEM CLS SZ 4 RIGHT
150401205	ATT CR FEM CLS SZ 5 RIGHT
150401206	ATT CR FEM CLS SZ 6 RIGHT
150401207	ATT CR FEM CLS SZ 7 RIGHT
150401208	ATT CR FEM CLS SZ 8 RIGHT
150401223	ATT CR FEM CLS SZ 3 N RIGHT
150401224	ATT CR FEM CLS SZ 4 N RIGHT
150401225	ATT CR FEM CLS SZ 5 N RIGHT
150401226	ATT CR FEM CLS SZ 6 N RIGHT
150401101	ATT CR FEM CLS SZ 1 LEFT

150401102	ATT CR FEM CLS SZ 2 LEFT
150401201	ATT CR FEM CLS SZ 1 RIGHT
150401202	ATT CR FEM CLS SZ 2 RIGHT
150401109	ATT CR FEM CLS SZ 9 LEFT
150401110	ATT CR FEM CLS SZ 10 LEFT
150401209	ATT CR FEM CLS SZ 9 RIGHT
150401210	ATT CR FEM CLS SZ 10 RIGHT
150600001	ATTUNE FB TIB BASE SZ 1 CEM
150600002	ATTUNE FB TIB BASE SZ 2 CEM
150600003	ATTUNE FB TIB BASE SZ 3 CEM
150600004	ATTUNE FB TIB BASE SZ 4 CEM
150600005	ATTUNE FB TIB BASE SZ 5 CEM
150600006	ATTUNE FB TIB BASE SZ 6 CEM
150600007	ATTUNE FB TIB BASE SZ 7 CEM
150600008	ATTUNE FB TIB BASE SZ 8 CEM
150600009	ATTUNE FB TIB BASE SZ 9 CEM
150600010	ATTUNE FB TIB BASE SZ 10 CEM
150610001	ATTUNE RP TIB BASE SZ 1 CEM
150610002	ATTUNE RP TIB BASE SZ 2 CEM
150610003	ATTUNE RP TIB BASE SZ 3 CEM
150610004	ATTUNE RP TIB BASE SZ 4 CEM
150610005	ATTUNE RP TIB BASE SZ 5 CEM
150610006	ATTUNE RP TIB BASE SZ 6 CEM
150610007	ATTUNE RP TIB BASE SZ 7 CEM
150610008	ATTUNE RP TIB BASE SZ 8 CEM
150610009	ATTUNE RP TIB BASE SZ 9 CEM
150610010	ATTUNE RP TIB BASE SZ 10 CEM
151620105	ATTUNE CR FB INSRT SZ 1 5 mm
151620106	ATTUNE CR FB INSRT SZ 1 6 mm
151620107	ATTUNE CR FB INSRT SZ 1 7 mm
151620108	ATTUNE CR FB INSRT SZ 1 8 mm
151620110	ATTUNE CR FB INSRT SZ 1 10 mm
151620112	ATTUNE CR FB INSRT SZ 1 12 mm
151620114	ATTUNE CR FB INSRT SZ 1 14 mm
151620116	ATTUNE CR FB INSRT SZ 1 16 mm
151620205	ATTUNE CR FB INSRT SZ 2 5 mm
151620206	ATTUNE CR FB INSRT SZ 2 6 mm
151620207	ATTUNE CR FB INSRT SZ 2 7 mm
151620208	ATTUNE CR FB INSRT SZ 2 8 mm
151620210	ATTUNE CR FB INSRT SZ 2 10 mm
151620212	ATTUNE CR FB INSRT SZ 2 12 mm
151620214	ATTUNE CR FB INSRT SZ 2 14 mm
151620216	ATTUNE CR FB INSRT SZ 2 16 mm

151620305	ATTUNE CR FB INSRT SZ 3 5MM
151620306	ATTUNE CR FB INSRT SZ 3 6MM
151620307	ATTUNE CR FB INSRT SZ 3 7MM
151620308	ATTUNE CR FB INSRT SZ 3 8MM
151620310	ATTUNE CR FB INSRT SZ 3 10MM
151620312	ATTUNE CR FB INSRT SZ 3 12MM
151620314	ATTUNE CR FB INSRT SZ 3 14MM
151620316	ATTUNE CR FB INSRT SZ 3 16MM
151620405	ATTUNE CR FB INSRT SZ 4 5MM
151620406	ATTUNE CR FB INSRT SZ 4 6MM
151620407	ATTUNE CR FB INSRT SZ 4 7MM
151620408	ATTUNE CR FB INSRT SZ 4 8MM
151620410	ATTUNE CR FB INSRT SZ 4 10MM
151620412	ATTUNE CR FB INSRT SZ 4 12MM
151620414	ATTUNE CR FB INSRT SZ 4 14MM
151620416	ATTUNE CR FB INSRT SZ 4 16MM
151620505	ATTUNE CR FB INSRT SZ 5 5MM
151620506	ATTUNE CR FB INSRT SZ 5 6MM
151620507	ATTUNE CR FB INSRT SZ 5 7MM
151620508	ATTUNE CR FB INSRT SZ 5 8MM
151620510	ATTUNE CR FB INSRT SZ 5 10MM
151620512	ATTUNE CR FB INSRT SZ 5 12MM
151620514	ATTUNE CR FB INSRT SZ 5 14MM
151620516	ATTUNE CR FB INSRT SZ 5 16MM
151620605	ATTUNE CR FB INSRT SZ 6 5MM
151620606	ATTUNE CR FB INSRT SZ 6 6MM
151620607	ATTUNE CR FB INSRT SZ 6 7MM
151620608	ATTUNE CR FB INSRT SZ 6 8MM
151620610	ATTUNE CR FB INSRT SZ 6 10MM
151620612	ATTUNE CR FB INSRT SZ 6 12MM
151620614	ATTUNE CR FB INSRT SZ 6 14MM
151620616	ATTUNE CR FB INSRT SZ 6 16MM
151620705	ATTUNE CR FB INSRT SZ 7 5MM
151620706	ATTUNE CR FB INSRT SZ 7 6MM
151620707	ATTUNE CR FB INSRT SZ 7 7MM
151620708	ATTUNE CR FB INSRT SZ 7 8MM
151620710	ATTUNE CR FB INSRT SZ 7 10MM
151620712	ATTUNE CR FB INSRT SZ 7 12MM
151620714	ATTUNE CR FB INSRT SZ 7 14MM
151620716	ATTUNE CR FB INSRT SZ 7 16MM
151620805	ATTUNE CR FB INSRT SZ 8 5MM
151620806	ATTUNE CR FB INSRT SZ 8 6MM
151620807	ATTUNE CR FB INSRT SZ 8 7MM

151620808	ATTUNE CR FB INSRT SZ 8 8MM
151620810	ATTUNE CR FB INSRT SZ 8 10MM
151620812	ATTUNE CR FB INSRT SZ 8 12MM
151620814	ATTUNE CR FB INSRT SZ 8 14MM
151620816	ATTUNE CR FB INSRT SZ 8 16MM
151620905	ATTUNE CR FB INSRT SZ 9 5 mm
151620906	ATTUNE CR FB INSRT SZ 9 6 mm
151620907	ATTUNE CR FB INSRT SZ 9 7 mm
151620908	ATTUNE CR FB INSRT SZ 9 8 mm
151620910	ATTUNE CR FB INSRT SZ 9 10 mm
151620912	ATTUNE CR FB INSRT SZ 9 12 mm
151620914	ATTUNE CR FB INSRT SZ 9 14 mm
151620916	ATTUNE CR FB INSRT SZ 9 16 mm
151621005	ATTUNE CR FB INSRT SZ 10 5 m
151621006	ATTUNE CR FB INSRT SZ 10 6 m
151621007	ATTUNE CR FB INSRT SZ 10 7 m
151621008	ATTUNE CR FB INSRT SZ 10 8 m
151621010	ATTUNE CR FB INSRT SZ 10 10 m
151621012	ATTUNE CR FB INSRT SZ 10 12 m
151621014	ATTUNE CR FB INSRT SZ 10 14 m
151621016	ATTUNE CR FB INSRT SZ 10 16 m
151640105	ATTUNE PS FB INSRT SZ 1 5 mm
151640106	ATTUNE PS FB INSRT SZ 1 6 mm
151640107	ATTUNE PS FB INSRT SZ 1 7 mm
151640108	ATTUNE PS FB INSRT SZ 1 8 mm
151640110	ATTUNE PS FB INSRT SZ 1 10 mm
151640112	ATTUNE PS FB INSRT SZ 1 12 mm
151640114	ATTUNE PS FB INSRT SZ 1 14 mm
151640116	ATTUNE PS FB INSRT SZ 1 16 mm
151640118	ATTUNE PS FB INSRT SZ 1 18 mm
151640205	ATTUNE PS FB INSRT SZ 2 5 mm
151640206	ATTUNE PS FB INSRT SZ 2 6 mm
151640207	ATTUNE PS FB INSRT SZ 2 7 mm
151640208	ATTUNE PS FB INSRT SZ 2 8 mm
151640210	ATTUNE PS FB INSRT SZ 2 10 mm
151640212	ATTUNE PS FB INSRT SZ 2 12 mm
151640214	ATTUNE PS FB INSRT SZ 2 14 mm
151640216	ATTUNE PS FB INSRT SZ 2 16 mm
151640218	ATTUNE PS FB INSRT SZ 2 18 mm
151640305	ATTUNE PS FB INSRT SZ 3 5MM
151640306	ATTUNE PS FB INSRT SZ 3 6MM
151640307	ATTUNE PS FB INSRT SZ 3 7MM
151640308	ATTUNE PS FB INSRT SZ 3 8MM

151640310	ATTUNE PS FB INSRT SZ 3 10MM
151640312	ATTUNE PS FB INSRT SZ 3 12MM
151640314	ATTUNE PS FB INSRT SZ 3 14MM
151640316	ATTUNE PS FB INSRT SZ 3 16MM
151640318	ATTUNE PS FB INSRT SZ 3 18MM
151640320	ATTUNE PS FB INSRT SZ 3 20MM
151640405	ATTUNE PS FB INSRT SZ 4 5MM
151640406	ATTUNE PS FB INSRT SZ 4 6MM
151640407	ATTUNE PS FB INSRT SZ 4 7MM
151640408	ATTUNE PS FB INSRT SZ 4 8MM
151640410	ATTUNE PS FB INSRT SZ 4 10MM
151640412	ATTUNE PS FB INSRT SZ 4 12MM
151640414	ATTUNE PS FB INSRT SZ 4 14MM
151640416	ATTUNE PS FB INSRT SZ 4 16MM
151640418	ATTUNE PS FB INSRT SZ 4 18MM
151640420	ATTUNE PS FB INSRT SZ 4 20MM
151640505	ATTUNE PS FB INSRT SZ 5 5MM
151640506	ATTUNE PS FB INSRT SZ 5 6MM
151640507	ATTUNE PS FB INSRT SZ 5 7MM
151640508	ATTUNE PS FB INSRT SZ 5 8MM
151640510	ATTUNE PS FB INSRT SZ 5 10MM
151640512	ATTUNE PS FB INSRT SZ 5 12MM
151640514	ATTUNE PS FB INSRT SZ 5 14MM
151640516	ATTUNE PS FB INSRT SZ 5 16MM
151640518	ATTUNE PS FB INSRT SZ 5 18MM
151640520	ATTUNE PS FB INSRT SZ 5 20MM
151640605	ATTUNE PS FB INSRT SZ 6 5MM
151640606	ATTUNE PS FB INSRT SZ 6 6MM
151640607	ATTUNE PS FB INSRT SZ 6 7MM
151640608	ATTUNE PS FB INSRT SZ 6 8MM
151640610	ATTUNE PS FB INSRT SZ 6 10MM
151640612	ATTUNE PS FB INSRT SZ 6 12MM
151640614	ATTUNE PS FB INSRT SZ 6 14MM
151640616	ATTUNE PS FB INSRT SZ 6 16MM
151640618	ATTUNE PS FB INSRT SZ 6 18MM
151640620	ATTUNE PS FB INSRT SZ 6 20MM
151640705	ATTUNE PS FB INSRT SZ 7 5MM
151640706	ATTUNE PS FB INSRT SZ 7 6MM
151640707	ATTUNE PS FB INSRT SZ 7 7MM
151640708	ATTUNE PS FB INSRT SZ 7 8MM
151640710	ATTUNE PS FB INSRT SZ 7 10MM
151640712	ATTUNE PS FB INSRT SZ 7 12MM
151640714	ATTUNE PS FB INSRT SZ 7 14MM

151640716	ATTUNE PS FB INSRT SZ 7 16MM
151640718	ATTUNE PS FB INSRT SZ 7 18MM
151640720	ATTUNE PS FB INSRT SZ 7 20MM
151640805	ATTUNE PS FB INSRT SZ 8 5MM
151640806	ATTUNE PS FB INSRT SZ 8 6MM
151640807	ATTUNE PS FB INSRT SZ 8 7MM
151640808	ATTUNE PS FB INSRT SZ 8 8MM
151640810	ATTUNE PS FB INSRT SZ 8 10MM
151640812	ATTUNE PS FB INSRT SZ 8 12MM
151640814	ATTUNE PS FB INSRT SZ 8 14MM
151640816	ATTUNE PS FB INSRT SZ 8 16MM
151640818	ATTUNE PS FB INSRT SZ 8 18MM
151640820	ATTUNE PS FB INSRT SZ 8 20MM
151640905	ATTUNE PS FB INSRT SZ 9 5 mm
151640906	ATTUNE PS FB INSRT SZ 9 6 mm
151640907	ATTUNE PS FB INSRT SZ 9 7 mm
151640908	ATTUNE PS FB INSRT SZ 9 8 mm
151640910	ATTUNE PS FB INSRT SZ 9 10 mm
151640912	ATTUNE PS FB INSRT SZ 9 12 mm
151640914	ATTUNE PS FB INSRT SZ 9 14 mm
151640916	ATTUNE PS FB INSRT SZ 9 16 mm
151640918	ATTUNE PS FB INSRT SZ 9 18 mm
151640920	ATTUNE PS FB INSRT SZ 9 20 mm
151641005	ATTUNE PS FB INSRT SZ 10 5 mm
151641006	ATTUNE PS FB INSRT SZ 10 6 mm
151641007	ATTUNE PS FB INSRT SZ 10 7 mm
151641008	ATTUNE PS FB INSRT SZ 10 8 mm
151641010	ATTUNE PS FB INSRT SZ 10 10 mm
151641012	ATTUNE PS FB INSRT SZ 10 12 mm
151641014	ATTUNE PS FB INSRT SZ 10 14 mm
151641016	ATTUNE PS FB INSRT SZ 10 16 mm
151641018	ATTUNE PS FB INSRT SZ 10 18 mm
151641020	ATTUNE PS FB INSRT SZ 10 20 mm
151650105	ATTUNE PS RP INSRT SZ 1 5 mm
151650106	ATTUNE PS RP INSRT SZ 1 6 mm
151650107	ATTUNE PS RP INSRT SZ 1 7 mm
151650108	ATTUNE PS RP INSRT SZ 1 8 mm
151650110	ATTUNE PS RP INSRT SZ 1 10 mm
151650112	ATTUNE PS RP INSRT SZ 1 12 mm
151650114	ATTUNE PS RP INSRT SZ 1 14 mm
151650116	ATTUNE PS RP INSRT SZ 1 16 mm
151650118	ATTUNE PS RP INSRT SZ 1 18 mm
151650205	ATTUNE PS RP INSRT SZ 2 5 mm

151650206	ATTUNE PS RP INSRT SZ 2 6 mm
151650207	ATTUNE PS RP INSRT SZ 2 7 mm
151650208	ATTUNE PS RP INSRT SZ 2 8 mm
151650210	ATTUNE PS RP INSRT SZ 2 10 mm
151650212	ATTUNE PS RP INSRT SZ 2 12 mm
151650214	ATTUNE PS RP INSRT SZ 2 14 mm
151650216	ATTUNE PS RP INSRT SZ 2 16 mm
151650218	ATTUNE PS RP INSRT SZ 2 18 mm
151650305	ATTUNE PS RP INSRT SZ 3 5MM
151650306	ATTUNE PS RP INSRT SZ 3 6MM
151650307	ATTUNE PS RP INSRT SZ 3 7MM
151650308	ATTUNE PS RP INSRT SZ 3 8MM
151650310	ATTUNE PS RP INSRT SZ 3 10MM
151650312	ATTUNE PS RP INSRT SZ 3 12MM
151650314	ATTUNE PS RP INSRT SZ 3 14MM
151650316	ATTUNE PS RP INSRT SZ 3 16MM
151650318	ATTUNE PS RP INSRT SZ 3 18MM
151650320	ATTUNE PS RP INSRT SZ 3 20MM
151650405	ATTUNE PS RP INSRT SZ 4 5MM
151650406	ATTUNE PS RP INSRT SZ 4 6MM
151650407	ATTUNE PS RP INSRT SZ 4 7MM
151650408	ATTUNE PS RP INSRT SZ 4 8MM
151650410	ATTUNE PS RP INSRT SZ 4 10MM
151650412	ATTUNE PS RP INSRT SZ 4 12MM
151650414	ATTUNE PS RP INSRT SZ 4 14MM
151650416	ATTUNE PS RP INSRT SZ 4 16MM
151650418	ATTUNE PS RP INSRT SZ 4 18MM
151650420	ATTUNE PS RP INSRT SZ 4 20MM
151650505	ATTUNE PS RP INSRT SZ 5 5MM
151650506	ATTUNE PS RP INSRT SZ 5 6MM
151650507	ATTUNE PS RP INSRT SZ 5 7MM
151650508	ATTUNE PS RP INSRT SZ 5 8MM
151650510	ATTUNE PS RP INSRT SZ 5 10MM
151650512	ATTUNE PS RP INSRT SZ 5 12MM
151650514	ATTUNE PS RP INSRT SZ 5 14MM
151650516	ATTUNE PS RP INSRT SZ 5 16MM
151650518	ATTUNE PS RP INSRT SZ 5 18MM
151650520	ATTUNE PS RP INSRT SZ 5 20MM
151650605	ATTUNE PS RP INSRT SZ 6 5MM
151650606	ATTUNE PS RP INSRT SZ 6 6MM
151650607	ATTUNE PS RP INSRT SZ 6 7MM
151650608	ATTUNE PS RP INSRT SZ 6 8MM
151650610	ATTUNE PS RP INSRT SZ 6 10MM

151650612	ATTUNE PS RP INSRT SZ 6 12MM
151650614	ATTUNE PS RP INSRT SZ 6 14MM
151650616	ATTUNE PS RP INSRT SZ 6 16MM
151650618	ATTUNE PS RP INSRT SZ 6 18MM
151650620	ATTUNE PS RP INSRT SZ 6 20MM
151650705	ATTUNE PS RP INSRT SZ 7 5MM
151650706	ATTUNE PS RP INSRT SZ 7 6MM
151650707	ATTUNE PS RP INSRT SZ 7 7MM
151650708	ATTUNE PS RP INSRT SZ 7 8MM
151650710	ATTUNE PS RP INSRT SZ 7 10MM
151650712	ATTUNE PS RP INSRT SZ 7 12MM
151650714	ATTUNE PS RP INSRT SZ 7 14MM
151650716	ATTUNE PS RP INSRT SZ 7 16MM
151650718	ATTUNE PS RP INSRT SZ 7 18MM
151650720	ATTUNE PS RP INSRT SZ 7 20MM
151650805	ATTUNE PS RP INSRT SZ 8 5MM
151650806	ATTUNE PS RP INSRT SZ 8 6MM
151650807	ATTUNE PS RP INSRT SZ 8 7MM
151650808	ATTUNE PS RP INSRT SZ 8 8MM
151650810	ATTUNE PS RP INSRT SZ 8 10MM
151650812	ATTUNE PS RP INSRT SZ 8 12MM
151650814	ATTUNE PS RP INSRT SZ 8 14MM
151650816	ATTUNE PS RP INSRT SZ 8 16MM
151650818	ATTUNE PS RP INSRT SZ 8 18MM
151650820	ATTUNE PS RP INSRT SZ 8 20MM
151650905	ATTUNE PS RP INSRT SZ 9 5 mm
151650906	ATTUNE PS RP INSRT SZ 9 6 mm
151650907	ATTUNE PS RP INSRT SZ 9 7 mm
151650908	ATTUNE PS RP INSRT SZ 9 8 mm
151650910	ATTUNE PS RP INSRT SZ 9 10 mm
151650912	ATTUNE PS RP INSRT SZ 9 12 mm
151650914	ATTUNE PS RP INSRT SZ 9 14 mm
151650916	ATTUNE PS RP INSRT SZ 9 16 mm
151650918	ATTUNE PS RP INSRT SZ 9 18 mm
151650920	ATTUNE PS RP INSRT SZ 9 20 mm
151651005	ATTUNE PS RP INSRT SZ 10 5 mm
151651006	ATTUNE PS RP INSRT SZ 10 6 mm
151651007	ATTUNE PS RP INSRT SZ 10 7 mm
151651008	ATTUNE PS RP INSRT SZ 10 8 mm
151651010	ATTUNE PS RP INSRT SZ 10 10 mm
151651012	ATTUNE PS RP INSRT SZ 10 12 mm
151651014	ATTUNE PS RP INSRT SZ 10 14 mm
151651016	ATTUNE PS RP INSRT SZ 10 16 mm

151651018	ATTUNE PS RP INSRT SZ 10 18 mm
151651020	ATTUNE PS RP INSRT SZ 10 20 mm
151630105	ATTUNE CR RP INSRT SZ 1 5 mm
151630106	ATTUNE CR RP INSRT SZ 1 6 mm
151630107	ATTUNE CR RP INSRT SZ 1 7 mm
151630108	ATTUNE CR RP INSRT SZ 1 8 mm
151630110	ATTUNE CR RP INSRT SZ 1 10 mm
151630112	ATTUNE CR RP INSRT SZ 1 12 mm
151630114	ATTUNE CR RP INSRT SZ 1 14 mm
151630116	ATTUNE CR RP INSRT SZ 1 16 mm
151630205	ATTUNE CR RP INSRT SZ 2 5 mm
151630206	ATTUNE CR RP INSRT SZ 2 6 mm
151630207	ATTUNE CR RP INSRT SZ 2 7 mm
151630208	ATTUNE CR RP INSRT SZ 2 8 mm
151630210	ATTUNE CR RP INSRT SZ 2 10 mm
151630212	ATTUNE CR RP INSRT SZ 2 12 mm
151630214	ATTUNE CR RP INSRT SZ 2 14 mm
151630216	ATTUNE CR RP INSRT SZ 2 16 mm
151630305	ATTUNE CR RP INSRT SZ 3 5MM
151630306	ATTUNE CR RP INSRT SZ 3 6MM
151630307	ATTUNE CR RP INSRT SZ 3 7MM
151630308	ATTUNE CR RP INSRT SZ 3 8MM
151630310	ATTUNE CR RP INSRT SZ 3 10MM
151630312	ATTUNE CR RP INSRT SZ 3 12MM
151630314	ATTUNE CR RP INSRT SZ 3 14MM
151630316	ATTUNE CR RP INSRT SZ 3 16MM
151630405	ATTUNE CR RP INSRT SZ 4 5MM
151630406	ATTUNE CR RP INSRT SZ 4 6MM
151630407	ATTUNE CR RP INSRT SZ 4 7MM
151630408	ATTUNE CR RP INSRT SZ 4 8MM
151630410	ATTUNE CR RP INSRT SZ 4 10MM
151630412	ATTUNE CR RP INSRT SZ 4 12MM
151630414	ATTUNE CR RP INSRT SZ 4 14MM
151630416	ATTUNE CR RP INSRT SZ 4 16MM
151630505	ATTUNE CR RP INSRT SZ 5 5MM
151630506	ATTUNE CR RP INSRT SZ 5 6MM
151630507	ATTUNE CR RP INSRT SZ 5 7MM
151630508	ATTUNE CR RP INSRT SZ 5 8MM
151630510	ATTUNE CR RP INSRT SZ 5 10MM
151630512	ATTUNE CR RP INSRT SZ 5 12MM
151630514	ATTUNE CR RP INSRT SZ 5 14MM
151630516	ATTUNE CR RP INSRT SZ 5 16MM
151630605	ATTUNE CR RP INSRT SZ 6 5MM

151630606	ATTUNE CR RP INSRT SZ 6 6MM
151630607	ATTUNE CR RP INSRT SZ 6 7MM
151630608	ATTUNE CR RP INSRT SZ 6 8MM
151630610	ATTUNE CR RP INSRT SZ 6 10MM
151630612	ATTUNE CR RP INSRT SZ 6 12MM
151630614	ATTUNE CR RP INSRT SZ 6 14MM
151630616	ATTUNE CR RP INSRT SZ 6 16MM
151630705	ATTUNE CR RP INSRT SZ 7 5MM
151630706	ATTUNE CR RP INSRT SZ 7 6MM
151630707	ATTUNE CR RP INSRT SZ 7 7MM
151630708	ATTUNE CR RP INSRT SZ 7 8MM
151630710	ATTUNE CR RP INSRT SZ 7 10MM
151630712	ATTUNE CR RP INSRT SZ 7 12MM
151630714	ATTUNE CR RP INSRT SZ 7 14MM
151630716	ATTUNE CR RP INSRT SZ 7 16MM
151630805	ATTUNE CR RP INSRT SZ 8 5MM
151630806	ATTUNE CR RP INSRT SZ 8 6MM
151630807	ATTUNE CR RP INSRT SZ 8 7MM
151630808	ATTUNE CR RP INSRT SZ 8 8MM
151630810	ATTUNE CR RP INSRT SZ 8 10MM
151630812	ATTUNE CR RP INSRT SZ 8 12MM
151630814	ATTUNE CR RP INSRT SZ 8 14MM
151630816	ATTUNE CR RP INSRT SZ 8 16MM
151630905	ATTUNE CR RP INSRT SZ 9 5 mm
151630906	ATTUNE CR RP INSRT SZ 9 6 mm
151630907	ATTUNE CR RP INSRT SZ 9 7 mm
151630908	ATTUNE CR RP INSRT SZ 9 8 mm
151630910	ATTUNE CR RP INSRT SZ 9 10 mm
151630912	ATTUNE CR RP INSRT SZ 9 12 mm
151630914	ATTUNE CR RP INSRT SZ 9 14 mm
151630916	ATTUNE CR RP INSRT SZ 9 16 mm
151631005	ATTUNE CR RP INSRT SZ 10 5 m
151631006	ATTUNE CR RP INSRT SZ 10 6 m
151631007	ATTUNE CR RP INSRT SZ 10 7 m
151631008	ATTUNE CR RP INSRT SZ 10 8 m
151631010	ATTUNE CR RP INSRT SZ 10 10 m
151631012	ATTUNE CR RP INSRT SZ 10 12 m
151631014	ATTUNE CR RP INSRT SZ 10 14 m
151631016	ATTUNE CR RP INSRT SZ 10 16 m
151820029	ATTUNE Medialized Dome Patella 29 mm
151820032	ATTUNE Medialized Dome Patella 32 mm
151820035	ATTUNE Medialized Dome Patella 35 mm
151820038	ATTUNE Medialized Dome Patella 38 mm

151820041	ATTUNE Medialized Dome Patella 41 mm
102407100	Sigma Uni Femoral component Left Medial / Right Lateral Size 1
102407200	Sigma Uni Femoral component Left Medial / Right Lateral Size 2
102407300	Sigma Uni Femoral component Left Medial / Right Lateral Size 3
102407400	Sigma Uni Femoral component Left Medial / Right Lateral Size 4
102407500	Sigma Uni Femoral component Left Medial / Right Lateral Size 5
102407600	Sigma Uni Femoral component Left Medial / Right Lateral Size 6
102408100	Sigma UNI femoral component Right Medial / Left Lateral Size 1
102408200	Sigma UNI femoral component Right Medial / Left Lateral Size 2
102408300	Sigma UNI femoral component Right Medial / Left Lateral Size 3
102408400	Sigma UNI femoral component Right Medial / Left Lateral Size 4
102408500	Sigma UNI femoral component Right Medial / Left Lateral Size 5
102408600	Sigma UNI femoral component Right Medial / Left Lateral Size 6
102451100	Sigma Uni Tibial component FB Left Medial / Right Lateral Size 1
102451200	Sigma Uni Tibial component FB Left Medial / Right Lateral Size 2
102451300	Sigma Uni Tibial component FB Left Medial / Right Lateral Size 3
102451400	Sigma Uni Tibial component FB Left Medial / Right Lateral Size 4
102451500	Sigma Uni Tibial component FB Left Medial / Right Lateral Size 5
102451600	Sigma Uni Tibial component FB Left Medial / Right Lateral Size 6
102452100	Sigma UNI Tibial component FB Right Medial / Left Lateral Size 1
102452200	Sigma UNI Tibial component FB Right Medial / Left Lateral Size 2
102452300	Sigma UNI Tibial component FB Right Medial / Left Lateral Size 3
102452400	Sigma UNI Tibial component FB Right Medial / Left Lateral Size 4
102452500	Sigma UNI Tibial component FB Right Medial / Left Lateral Size 5
102452600	Sigma UNI Tibial component FB Right Medial / Left Lateral Size 6
102453107	Sigma UNI FB inserts Left Medial / Right Lateral Size 1 7
102453108	Sigma UNI FB inserts Left Medial / Right Lateral Size 1 8
102453109	Sigma UNI FB inserts Left Medial / Right Lateral Size 1 9
102453110	Sigma UNI FB inserts Left Medial / Right Lateral Size 1 10
102453111	Sigma UNI FB inserts Left Medial / Right Lateral Size 1 11
102453207	Sigma Uni FB inserts Left Medial / Right Lateral Size 2 7
102453208	Sigma Uni FB inserts Left Medial / Right Lateral Size 2 8
102453209	Sigma Uni FB inserts Left Medial / Right Lateral Size 2 9

102453210	Sigma Uni FB inserts Left Medial / Right Lateral Size 2 10
102453211	Sigma Uni FB inserts Left Medial / Right Lateral Size 2 11
102453307	Sigma UNI FB inserts Left Medial / Right Lateral Size 3 7
102453308	Sigma UNI FB inserts Left Medial / Right Lateral Size 3 8
102453309	Sigma UNI FB inserts Left Medial / Right Lateral Size 3 9
102453310	Sigma UNI FB inserts Left Medial / Right Lateral Size 3 10
102453311	Sigma UNI FB inserts Left Medial / Right Lateral Size 3 11
102453407	Sigma UNI FB inserts Left Medial / Right Lateral Size 4 7
102453408	Sigma UNI FB inserts Left Medial / Right Lateral Size 4 8
102453409	Sigma UNI FB inserts Left Medial / Right Lateral Size 4 9
102453410	Sigma UNI FB inserts Left Medial / Right Lateral Size 4 10
102453411	Sigma UNI FB inserts Left Medial / Right Lateral Size 4 11
102453507	Sigma UNI FB inserts Left Medial / Right Lateral Size 5 7
102453508	Sigma UNI FB inserts Left Medial / Right Lateral Size 5 8
102453509	Sigma UNI FB inserts Left Medial / Right Lateral Size 5 9
102453510	Sigma UNI FB inserts Left Medial / Right Lateral Size 5 10
102453511	Sigma UNI FB inserts Left Medial / Right Lateral Size 5 11
102453607	Sigma UNI FB Left Medial / Right Lateral Size 6 7
102453608	Sigma UNI FB Left Medial / Right Lateral Size 6 8
102453609	Sigma UNI FB Left Medial / Right Lateral Size 6 9
102453610	Sigma UNI FB Left Medial / Right Lateral Size 6 10
102453611	Sigma UNI FB Left Medial / Right Lateral Size 6 11
102454107	Sigma UNI FB Right Medial / Left Lateral Size 1 7
102454108	Sigma UNI FB Right Medial / Left Lateral Size 1 8
102454109	Sigma UNI FB Right Medial / Left Lateral Size 1 9
102454110	Sigma UNI FB Right Medial / Left Lateral Size 1 10
102454111	Sigma UNI FB Right Medial / Left Lateral Size 1 11
102454207	Sigma UNI FB Right Medial / Left Lateral Size 2 7
102454208	Sigma UNI FB Right Medial / Left Lateral Size 2 8
102454209	Sigma UNI FB Right Medial / Left Lateral Size 2 9
102454210	Sigma UNI FB Right Medial / Left Lateral Size 2 10
102454211	Sigma UNI FB Right Medial / Left Lateral Size 2 11
102454307	Sigma UNI FB inserts Right Medial / Left Lateral Size 3 7
102454308	Sigma UNI FB inserts Right Medial / Left Lateral Size 3 8
102454309	Sigma UNI FB inserts Right Medial / Left Lateral Size 3 9
102454310	Sigma UNI FB inserts Right Medial / Left Lateral Size 3 10

102454311	Sigma UNI FB inserts Right Medial / Left Lateral Size 3 11
102454407	Sigma UNI FB inserts Right Medial / Left Lateral Size 4 7
102454408	Sigma UNI FB inserts Right Medial / Left Lateral Size 4 8
102454409	Sigma UNI FB inserts Right Medial / Left Lateral Size 4 9
102454410	Sigma UNI FB inserts Right Medial / Left Lateral Size 4 10
102454411	Sigma UNI FB inserts Right Medial / Left Lateral Size 4 11
102454507	Sigma UNI FB inserts Right Medial / Left Lateral Size 5 7
102454508	Sigma UNI FB inserts Right Medial / Left Lateral Size 5 8
102454509	Sigma UNI FB inserts Right Medial / Left Lateral Size 5 9
102454510	Sigma UNI FB inserts Right Medial / Left Lateral Size 5 10
102454511	Sigma UNI FB inserts Right Medial / Left Lateral Size 5 11
102454607	Sigma UNI FB inserts Right Medial / Left Lateral Size 6 7
102454608	Sigma UNI FB inserts Right Medial / Left Lateral Size 6 8
102454609	Sigma UNI FB inserts Right Medial / Left Lateral Size 6 9
102454610	Sigma UNI FB inserts Right Medial / Left Lateral Size 6 10
102454611	Sigma UNI FB inserts Right Medial / Left Lateral Size 6 11