



SURGIVAL
g r u p o c o s i a s

Valencia 12th July 2004

By the present the company SURGIVAL CO S.A.U, informs that:

- ✓ It is possible to do a "Tomography" for the self-looking stems, bipolar heads and polyethylene cups
It does not exist any type of incompatibility between the tomographical technic (TAC) and the material used in the manufacturing of these products.
The only problem can be the resolution of the image due to the lack of clarity, unless that you are provided with the equipment of image processing.
- ✓ Regarding the "Magnetic Resonance" (NMR) there are not studies published which says that a Magnetic Resonance can damage, being perjudicial or cause secondary effects to the patient.
The only problems which can cause does not affect to the patient, but as it is written in our Postoperative Precautions:
 - There is a possibility of electromagnetic interference with certain rehabilitation systems based on the emission of high-frequency waves.
 - This product is detected by metal detection systems
 According to some articles published regarding this matter it is written that only on ferromagnetics materials as the stainless steel it is not convenient the use of this technic because it can gives problems on the image resolution.

Yours sincerely,

M^o José Hernández
Technical Responsible
SURGIVAL Co S.A.U.

SURGIVAL Co, S.A.U.
Parque Tecnológico
C/ Leonardo Da Vinci, 12 - 14
Teléfono. 96 136 62 66 - 96 131 80 50
Fax 96 131 80 95
46980 PATERNA (Valencia)

SURGIVAL Co, S.A.U. · Parque Tecnológico - C/ Leonardo Da Vinci, 12-14 · 46980 Paterna (Valencia) Spain
Tel. 00 34 96 131 80 50 / 96 136 62 66 · Fax 00 34 96 131 80 95 · E-mail: surgival@surgival.com · www.surgival.com